
HM-Fräser
TC-Cutter
Fraises en carbure

		Drehzahlen entsprechend der Materialhärten und Arbeitsteil grössen: Speeds according to the material hardnesses and working steep size: ↻_{opt.} 10'000 - 30'000 min⁻¹																				
		10	20	24	25	26	27	30	31	40	41	50	51	53	55	60	70	75	80	85	90	L10 L20 L55
Metall-alloys-alliages																						
Edelmetall-Legierungen semi-precious alloys alliages semi précieux	standard	fein fine fine	fein fine fine	superfein superfine superfine	grob coarse grosse	fein fine fine	superfein superfine superfine	superfein superfine superfine	fein fine fine	standard						superfein superfine superfine	standard					universell universal univer'sell
NEM-Legierungen non-precious alloys alliages non précieux	standard	fein fine fine	fein fine fine	superfein superfine superfine	grob coarse grosse	fein fine fine	superfein superfine superfine	superfein superfine superfine	fein fine fine	standard						superfein superfine superfine	standard					universell universal univer'sell
Modellguss model cast appliances squelettés	standard	fein fine fine		superfein superfine superfine									mittelgrob med.coarse moy.grosse				standard					
Titan titanium titane	standard		fein fine fine	superfein superfine superfine	grob coarse grosse	fein fine fine																
Gips-plaster-plâtre																						
Gips plaster plâtre	standard			superfein superfine superfine								grob coarse grosse	mittelgrob med.coarse moy.grosse		supergrob sup. coarse sup.grosse		standard		grob coarse grosse	supergrob sup. coarse sup.grosse	grob coarse grosse	universell universal univer'sell
Keramik-ceramics-céramique																						
Keramik ceramics céramique		fein fine fine						superfein superfine superfine	superfein superfine superfine													universell universal univer'sell
Kunststoff-resin-résin																						
Kunststoffe resins résine	standard	fein fine fine	fein fine fine	superfein superfine superfine								grob coarse grosse	mittelgrob med.coarse moy.grosse	feingrob fin. coarse fin.grosse	supergrob sup. coarse sup.grosse	superfein superfine superfine	standard		grob coarse grosse	supergrob sup. coarse sup.grosse		universell universal univer'sell
Weichbleibende Unterfütterung soft relinings rebasages moux																			Querrieb transverse transversale			
Composite composites composite								superfein superfine superfine	superfein superfine superfine													universell universal univer'sell
Maximal-Drehzahl Das Nichtbeachten der maximal zulässigen Drehzahl führt zu einem erhöhten Sicherheitsrisiko. Maximum-Speeds Non-adherence to the maximum permissible speeds increases the risk of accidents. Vitesse-Maximale Le dépassement de la vitesse de rotation maximale permise constitue un risque de sécurité élevé.																						
ISO Ø 1/10 mm																						
010 – 023	upm																					
025 – 045	5.000 – 40.000																					
050 – 080	15.000 – 30.000 15.000 – 20.000																					
Grobausarbeiten, Konturieren Rough trimming, contouring Façonnage, retouches																						
Feinausarbeiten, Glätten Fine elaboration, smoothing of surfaces Polissage des surfaces, corrections																						
Ausarbeiten, Konturieren Trimming, contouring Façonnage, retouches																						
Ausarbeiten, Konturieren Trimming, contouring Façonnage, retouches																						
Ausarbeiten, Konturieren Trimming, contouring Façonnage, retouches																						
Feinausarbeiten, Glätten Fine elaboration, smoothing of surfaces Polissage des surfaces, corrections																						
Feinausarbeiten, Glätten Fine elaboration, smoothing of surfaces Polissage des surfaces, corrections																						
Feinausarbeiten, Glätten Fine elaboration, smoothing of surfaces Polissage des surfaces, corrections																						
Verbundfördernde Strukturierung Active bonding structuring Préparer la surface																						
Verbundfördernde Strukturierung Active bonding structuring Préparer la surface																						
Grober Materialabtrag Bulk material reduction Rapide réduction																						
Grober Materialabtrag Bulk material reduction Rapide réduction																						
Grobausarbeiten, Konturieren Rough trimming, contouring Façonnage, retouches																						
Grober Materialabtrag Bulk material reduction Rapide réduction																						
Feinausarbeiten, Glätten Fine elaboration, smoothing of surfaces Polissage des surfaces, corrections																						
Grobausarbeiten, Konturieren Rough trimming, contouring Façonnage, retouches																						
Ausarbeiten Trimming Façonnage																						
Grober Materialabtrag Bulk material reduction Rapide réduction																						
Grober Materialabtrag Bulk material reduction Rapide réduction																						
Spannungsfreies Ausbetten Stress-free deflasking Pour l'élimination douce																						
Linkshändfräse For left hand use Fraise pour gaucher																						

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Please Note: Using cutters marked with the indentation codes 40 - 85 the grinding movements must be done under constant pressure operating the instrument at the recommended speed.
Danger of Injury – soft, timid or hesitant use of large intended instruments causes rock or toss vibration with the possibility of a shaft breakage.

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