## **zubler** Vacuum Suction Systems

Appropriate solutions for any laboratory





www.zubler.de

#### Grinding without suction?

Laboratory technicians would occasionally prefer to turn off their vacuum during grinding.

Even during milling or grinding jobs that scarcely involve any reduction of material, fine dust is always created. You cannot see them, and yet they are whirling and hovering across the room. With your head located close to the instrument tip, you will inevitably inhale particles in high concentrations even when performing grinding jobs that you may consider perfectly harmless. It is precisely these fine dust particles that will penetrate deeply into the lungs with ease. Common particles from alveolar bone (< 5  $\mu$ m in diameter) can become lodged in the lungs and turn into a source of chronic asthmatic conditions and severe lung disease.

Think about it! We are talking about your health. You perform these jobs almost daily.

## Occupational health and safety are non-negotiable

If dental laboratories fail to use suction, most dust that is generated will clearly exceed occupational exposure limits.

Suction to eliminate dusts is invariably required by German regulations governing the use of hazardous substances. Furthermore, all vacuum systems must be mandatorily tested for effectiveness at least once a year.

#### We support you with our experience

Dust exposure can be dramatically reduced by the use of optimized funnel systems and state-of-the-art vacuum equipment.

Vacuum systems are particularly effective if they convey exhaust air to the outside of the building. This will also improve the indoor climate throughout the laboratory.

Zubler offers a number of vacuum options, so you can always select a solution that is optimized for the needs of your laboratory.

#### Turn on your vacuum – always!

Even vacuum systems not using the latest technology will usually reduce dust levels around the workspace. However, you should make an effort to improve the situation before there is a problem. Single suction unit

Multi suction unit

**Central suction system** 

Accessories

## Dust in dental technology: a big occupational hazard

Dental technicians are especially at risk. According to statistics by the BGFE (German Occupational Society for Precision and Electrical Engineering), occupational skin and lung diseases are considerably more prevalent among dental technicians than the BGFE average. Fine airborne particles are especially dangerous. Highly detrimental effects on the skin and respiratory organs have also been identified for dental gypsum and investment materials. The same is true of non-classified particles with general exposure limits.



## Product range selection

#### Vac. ports

#### Single suction unit



#### **Multi suction unit**



## FZ1 VARIOmaster®

**V4000** 

V5000 V6000

V7000

Our single-position units are operated with filter bags. Various models with different filters and suction motors are offered, which are optimized for different dust

loads and applications.

FZ VARIO

All of our multiple-position units are equipped with non-brush motors and permanent filter systems. The dust is collected in a container. Although these systems are more expensive, they are highly economical in the long run.

#### **Central suction system**



#### Accessories

Suction hood
rge vacuum funne
Adjustable pipe
nall vacuum funnel
Gold filter insert
Rectangular pipe
Intake system

Central vacuum systems offer the best working environment. They are the cleanest and quietest performers. Whenever a solution for a larger number of vacuum positions needs to be found, it is a good idea to analyze whether a central system might be economical and could be accommodated.



## 2 (4)

>6

El. pov min.	<b>ver</b> (W) max.	Air volu min.	<b>ıme</b> (I/s) max.	Noise min.	(dB(A)) max.	Ű	Anaust of	utdoors traust in Fi	doors iterbag	ter carti	idge stoonbring	ush motor	notor of the pass cooling		
						_		_							
250	650	25	36	55	65	-		•	-	•	-	-			Ŀ,
300	1100	25	45	62	77	•	-	-		-	-	•			
250	650	25	38	52	61	-	-	-			-	-			io
150	700	20	40	54	66	-	-			-	-	-			uct
													Pag	e <b>6-9</b>	Single s
200	1200	16	65	48	67	•	•	-	•	-	•	•			Ś
420	2200	20	90	49	60	•	<u> </u>	-	•	-	•		Page 1	0-15	Multi suction unit
600	2500	20	300	(52)	(71)	•			•			•	Page 1	6-21	Central suction system
R11 R10 AV MB P10 FG ST0	500 600 1000 1000 1000 1000 12	Vacuum h Monomer Preparatio Monomer Polishing Base Vacuum p Cleaning a	bath n station box unit lug acces.										Page 2	2-26	Accessories

## V4000

#### **Mobile suction**

#### Workspace vacuum for installation and portable use

- Compact unit featuring an exchangeable filter door system.
- Continuous electronic control of vacuum power.
- Automatic on-off with connection to device.
- Switches off when the filter is full.
- Different bag filters and activated carbon filters are offered as accessories.





## Filter technology of the **V4000**

All **V4000** units feature a dual filter system. Dirt is collected in the filter bag, which allows for easy removal and clean disposal without any dust escaping. A fine filter cartridge offers extra protection. It has longer replacement intervals.

Filter quality complies with dust category M as defined in EN 60335-2-69.



#### Filters

A side from the standard filter bag, a large filter bag is available for increased efficiency and extended service life in high-dirt

applications. The large bag only works in conjunction with the wide **V4050** filter door, which is mutually compatible with the **V4000** unit to optimize versatility.



#### Activated charcoal filter

An activated charcoal filter can be inserted instead of the filter cartridge. Activated charcoal is capable of absorbing toxic vapours (for example, when resins are processed or monomers mixed). The activated charcoal filter is supplied as a refillable cartridge. Refill packs are available for convenient, environmentally friendly and inexpensive replacement of the charcoal.



#### V4011 built-in version

The built-in version is suspended under the benchtop. With the mounting plate and matching rail available for this purpose, the unit can be readily installed and removed. Even the mobile unit is easily convertible to a built-in version. Just replace the carrying handle with the mounting rail (optional accessory).



#### **Technical Data**

**V4000** 

Dimensions:		
Width:	small / wide door	265/ 375mm
Height:	handle / mounting rail	435/ 415mm
Depth:		400mm
Weight:	small / wide door	20/25kg
Power input:	suction unit	max. 700W
	autom. socket	max. 450W
Suction power:		18-35 l/s
Noise level:		55-65 dB(A)

## V5**000**

#### **Compact mobile suction**

#### Device/workspace vacuum. Offers a brushless motor.

- Compact unit featuring an exchangeable filter door system.
- Collectorless motor with continuous operation properties.
- Different bag filters and activated carbon filters are offered as accessories.



## V6000

#### Workbench suction

## Workspace vacuum. Offers low-noise operation.

- Narrow tower housing with a large filter volume.
- Front access for filter maintenance through a caddy bay.
- Very low noise by high-grade acoustic isolation.
- Different bag filters and activated charcoal filters are offered as accessories.





#### **Tower suction**

## **Device/workspace vacuum.** Offers excellent suction and prolonged operation.

- Narrow tower housing with a large filter volume.
- Front access for filter maintenance through a caddy bay.
- Brushless motor.
- Different bag filters and activated charcoal filters are offered as accessories.



#### **Application examples**



Tee





#### **Applications**

Our single suction units are highly versatile, as various models with different filter systems and motors are available. For low-dust applica-tions, they can also be used with manually controlled or automated milling devices (e.g. for zirconia objects). Appropriate suction units can also be selected for any workshop or laboratory applicati-ons involving dry dusts. Applications of this type may inclu-de various modes of processing minerals, resins or metals in diffe-rent sectors of industry.

<b>Technical Data</b>		V5000	V6000	V7000	
Dimensions:					
Width:		375mm	200mm	200mm	ر س
Height:	handle / mounting rail	440/ 420mm	650mm	650mm	rie
Depth:		450mm	600mm	600mm	SSO
Weight:		20kg	31kg	31kg	Ce Ce
Power input:	suction unit	max. 1100W	650W	700W	Ă
	autom. socket	- / / /	450W	500W	
Suction power:		25-45 l/s	25-38 l/s	20-40 l/s	
Noise level:		62-77 dB(A)	52-61 dB(A)	54-66 dB(A)	

## FZ1 VARIOmaster®

#### 2-Station-unit

#### Good reasons for an FZ1

- Low noise and low maintenance (thanks to a brushless vacuum motor offering bypass air cooling and sophisticated acoustic isolation of the exhaust stream).
- Permanent filter with automatic filter cleaning system. No filter bags.
- Separately exchangeable and inexpensive filter cartridges for particle separation and fine filtration.
- High-performance vacuum for use with belt sanders and sandblasters.
- Can be expanded via external position openers to accommodate up to four alternatively used vacuum positions.
- More room at the workspace.
- Comes ready with an exhaust connection, but recirculation is also an option.
- High-grade components and modular design for long service life and easy maintenance.



## Vibration cleaned filter cartridge

How well a vacuum system performs will depend on the quality of the filter material and on the efficiency of the cleaning function of the filter. In addition to the brushless vacuum motor, the **FZ1 VARIOmaster**<sup>®</sup> utilizes a highly sophisticated filter technology. The filters are arranged such that they can vibrate freely inside the housing. High-grade vibration equipment creates extremely powerful vibration, thus ensuring that dust layers are eliminated from the filter surface time after time.



FZ1 VARIOmaster®

FZ1 VARIOmaster® is suitable for universal use in laboratories and workshops. A highly flexible concept allows for connection of four vacuum positions, to be controlled through individual modules. Vacuum positions are supplied as requested by the customer. This robust and high-performance suction unit is equally well suited for workspaces and for high-dust laboratory equipment. Laboratories in dental offices, or small-sized dental laboratories including only one or two technicians, may be able to cover their needs for device or workspace suction with a single FZ1 VARIOmaster<sup>®</sup> unit. Particularly comfortable is its alternative mode of operation, which allows different vacuum power levels to be defined for each position connected. The VARIO control system is simultaneously active, maintaining a constant air volume at all times, regardless of the condition of the filter.

#### Automatic vacuum port openers

AP501 modules will allow the **FZ1** VARIOmaster<sup>®</sup> to be configured for up to 4 vacuum positions. Each AP501 offers connectivity for any universal controllers or laboratory equipment with up to 500W of power consumption. Its electronic control offers sensitivity adjustment to match the

connected (dust-generating) equipment, switches the vacuum position valve, and connects to the suction unit via a control cable.

Both the valve and the control unit are always located at the workspace. This makes installation very easy and transparent. Several vacuum positions can be connected to the same pipeline.



#### **Dust container** (No filter bags)

The high performance of this filter system comes into its own in very high-dirt areas. The use of a longlasting filter cartridge eliminates the need for filter bags, thus greatly reducing operating costs.

All dust coming from the connected vacuum positions is collected in a container, which is emptied as needed.



#### **Technical Data**

#### FZ1 VARIOmaster®

Dimensions: Width x Height x Depth:	224 x 760 x 600mm
Weight:	38kg
Power input:	max. 1100W
Motor:	brushless
Suction power:	approx. 60 l/s
Dust container volume:	5
Filter quality:	Dust class M
Filter surface:	1,5 m²

#### "VARIO"- the system with better motors

Brushless motors have a long ser-

The motor of the **FZ1 VARIOmaster**<sup>®</sup> is cooled by a ventilation circuit that bypasses (thus being separated from) the vacuum stream.

**Multi suction units** 

## FZ2 VARIOmatic<sup>®</sup>

#### 4-Station-unit

#### Good reasons for an FZ2

- Two vacuum turbines (three-phase AC) for low noise and low maintenance. No brushes.
- Permanent filter with automatic filter cleaning system. No filter bags.
- Pneumatic filter reverse flow for high-performance suction of major dust loads.
- Can be expanded via external position openers to accommodate up to six vacuum positions.
- More room at the workspace.
- Exhaust installation for a better indoor climate.





FZ2 VARIOmatic®

This unit was designed specifically for use in small dental laboratories. It offers a particularly reliable and clean vacuum solution even for high-dust preparatory workspaces and for functional equipment such as belt sanders. Office laboratories or small-sized regular laboratories may be able to cover all requirements for device or workspace suction by installing a single **FZ2 VARIOmatic**<sup>®</sup> as a mini-central.

When a handpiece or other connected tool is activated, each single vacuum position opens automatically (and regardless of any other positions connected). Once the tool is deactivated, it will close with a short delay to eliminate any residual dust (delayed-stop function).

**FZ2 VARIOmatic**<sup>®</sup> is continuously on standby, building up vacuum outputs as required without delay. Control of the vacuum turbines will automatically adjust to the number of open vacuum positions.



FZ2 VARIOmatic® Dentallab Engelhardt + Schnierer - Ulm/Germany 🔺

## Correct planning of vacuum installation

**FZ2 VARIOmatic**<sup>®</sup> will integrate with any furniture line. Matching front panels are mounted on request. Zubler will be happy to advise you and support your planning. Our experience will help you to find the correct vacuum solution for your laboratory.

Installation and maintenance will be greatly facilitated if the connections of the vacuum system are well accessible. Exhaust air should take the shortest path to the outside. Taking advantage of these options will greatly improve your lab environment and reduce noise even further.

#### **Expansion modules**

**FZ2 VARIOmatic**<sup>®</sup> can even be installed to cover 5 or 6 vacuum positions if AP500-D, AP501 or FP1000-D modules are additionally used. AP500-D and AP501 will connect to universal controllers, offering electronics with sensitivity adjustment. FP1000-D will connect to laboratory devices with up to 1000 W of power consumption. The exhaust air from the **FZ2 VARIOmatic**<sup>®</sup> must be conveyed to the outside of the building. The entire stream of exhaust air is blasted outside through a piping connection 75 mm in diameter, which can be linked up with commercially available HT pipes.

The optimal solution is the standard variation shown by the connections of the FZ2 VARIOmatic® to 4 hand nozzle working positions close to each other. All working positions can be easily and quickly connected with a Ø=40mm hose. The vacuum system is easily reached from the back side for maintenance through the free area at the base and the free area under the flower pot.

The **FZ2 VARIOmatic**<sup>®</sup> is especially suitable for use where a very high amount of dust is produced, e.g. in an instrument room. Another unit (i.e. a fifth vacuum position) can be connected externally through an electronic switchbox with an additional valve. Pipes rather than individual flexible tubes should be installed if longish distances need to be covered.

Installations of this type are slightly more complex. At the same time, they result in a more compact system that is cleaner and offers better flow efficiency. In this example, the only direct connection via a flexible tube was established for the nearby sandblaster.

Exhaust

The **FZ2 VARIOmatic**<sup>®</sup> can, in some cases, e.g. in Practices and small laboratories, cover the entire requirement for vacuum positions. If some working positions are alternatively used, the connection of up to 6 vacuum positions is possible. All of the dirt from the laboratory would be collected in one container. Blowing the exhaust to the outside gives an optimal solution concerning noise, air quality, room climate and economy.

This example shows another situation in which a shared pipeline for all vacuum positions would be appropriate.







#### Addressing problems at their roots

Collectorless motors have a long service life. Zubler offer a three-year warranty. Other advantages of these motors include:

- No need to replace brushes.
- No carbon dust in the exhaust air.
- Operation with less noise and vibration.



Implementing motors of this type is a complex task, but the returns become noticeable very quickly in vacuum systems used as a centre for several positions.

FZ2 VARIOmatic® is equipped with two brushless motors.



#### Air performance and power consumption: filter technology rules

How well a vacuum system performs will depend on the size and guality of the filter material and on the efficiency of the cleaning function of the filter. Much of the filter technology used in competing products leaves something to be desired. FZ2 VARIOmatic<sup>®</sup>, by contrast, has been supplied with a filter cartridge encompassing a very large filter area (2.3 m<sup>2</sup>). This feature will ensure a high degree of air permeability.

Dusts are retained in the fine mesh structure of the filter material. The system features a pneumatic reverse flow, which will ensure that this mesh is continuously blasted clean.



FZ2 VARIOmatic<sup>®</sup> Filter cartridge

#### **Technical Data**

FZ2 VARIOmatic®

Dimensions: Width x Height x	Depth:	334 x 810 x 600mm
Weight:		72kg
Power input:	suction unit:	1900W
	autom. sockets:	4x 450W
Motor:		brushless 2pcs.
Suction power:		approx. 80 l/s
Dust container volume:		12
Filter surface:		2,3m <sup>2</sup>
Necessary pressure:		min. 5bar

#### Effective filter equipment does have a footprint

Suction units should preferably be for high air performance will occu-py considerable space. High-per-formance blowers are either bulky, or they operate at high speed. In other words, smaller blowers are noisier. Good acoustic isolation can only work if certain surface and volume requirements are met.

We have delivered what is technounderstanding that vacuum equipfind an appropriate solution that is powerful and will keep you com-fortable in the long run. We put

#### **FZ VARIO** Central suction system



#### Good reasons for a "center"

#### Improved air quality

- Exhaust air is conveyed to the outside, which is very beneficial for the indoor climate.
- When resins are processed, all toxic and unpleasant vapours are eliminated, and the working environment is cleared of residual/microscopic dusts.
- Air change is promoted.

#### Less noise

• No vacuum motor around the workspace.

#### More space

 As no suction units and filter housings are located around the workspace, the technician has more space for tools/instruments and more freedom of movement.

#### Greater cleanliness

- No exchanging of filters in the laboratory, and dust is evacuated through the pipeline to the vacuum unit.
- Cleaning of the laboratory is made easier by preinstalled vacuum outlets.

#### Smaller overhead

• High-efficiency large blower motors for lower energy consumption.

The FZ VARIO central vacuum system with its accessory intake funnels is a solution without compromise. With its efficient vacuum motors and filtration systems, available in different output classes and combinations, the Zubler central vacuum systems are suitable for providing suction at all dedicated workspaces for handpiece or resin processing, or indeed for any dustproducing devices, in dental laboratories of any size. The advantages of a central vacuum system become particularly obvious when all dustgenerating devices and workspaces in the laboratory benefit from central suction. Devices with internal blower motors can be retrofitted for connectivity to the central vacuum system.

The VARIO system controls the air volume stream, continuously adapting the suction depending on the number of open vacuum positions. Compensation for filter pollution is also offered. The filter system will perform autocleaning at defined intervals, such that a high level of air permeability is sustained continuously. The motor seldom reaches its full speed and only does so briefly. These features will ensure a long service life and keep maintenance requirements very low. Laboratories using the **FZ VARIO** will observe all pertinent German regulations (German Federation of Commercial and Industrial Compensation Societies, Trade Supervisory Authority).



#### Automatic vacuum port openers

Automatic vacuum position openers are standard features of any Zubler central vacuum system. They contribute significantly to the economical operation of the system, since the vacuum will be provided to the workspace only while dust is actually being generated. When a handpiece or other connected tool is activated (e.g. a sandblaster), the vacuum position opens automatically. Once the tool is deactivated, it will close with a short delay to eliminate any residual dust (delayedstop function).

This process of automatically linking the dust-generating device to the electronic control of the vacuum position opener will ensure that your health is protected and that pertinent legal requirements for occupational safety and health are met (e.g. Germany).

# Multi suction units

Single suction unit

### Modules

#### AP 500

AP500 is a vacuum position opener for all common universal controllers. Its response sensitivity can be automatically calibrated or manually adjusted. Being suitable for any laboratory devices up to 500W, this unit is the most flexible and highgrade controller for workspaces needing suction.

#### FP-D and FP 1000

FP-D is a vacuum position opener not featuring a delayed-stop function. It is designed for use with laboratory devices offering a directconnection outlet.

FP1000 is a vacuum position featuring a delayed-stop function. It is designed for use with laboratory devices not accepting standby current and not featuring a vacuum outlet. This unit does not have suitable electronics for handpiece positions.

#### FP-M

FP-M is a vacuum position opener for manual activation. It can be used at workspaces where there is no electrical device capable of automatically opening the vacuum position, and which are exclusively designed for manual operation (e. g. mixing of monomer).

#### **Special modules**

#### Module F500 for Freuding

Special table inserts have been designed to match Freuding's F20 and F40 furniture lines. In addition to the electronic vacuum position opener, these inserts feature a pullout suction canal with an integrated sound damper. Zubler's intake funnels R1200 and R1250 can be used with this system.



#### Module K500 for KaVo

With Ergospace furniture by KaVo, instead of using the vacuum drawer, the Zubler module is simply advanced into the existing support frame.



Custom-tailored solutions are also available for Masterspace and Flexspace workstations. Zubler's membrane keyboard is worked into KaVo's original panel.



Precious-metal dust should be collected as efficiently as possible in order to recover as much valuable material. The cost of precious metal recycling can be minimized by ensuring a high precious-metal content of the recyclable material. The Zubler two-circuit technology uses a separate filtering system for all positions where the dust generated is likely to be (1) dust from a pure material such as dental stone, acrylic or steel or (2) dust containing residual precious metal with only a low proportion of other dust. This requires installation of two separate pipe systems, which can be connected to the respective filtration system.



#### Separation of precious metals

- High-performance suction and effective health protection during precious-metal jobs.
- Superior performance to filter inserts in the intake funnel (which become rapidly clogged)
- Centralized detection of precious-metal dusts (independent of the technician's care and attention).
- Larger overall filter area for better performance of the vacuum system.
- Direct disposal of bulk dirt.
- Reduced separation costs.

#### Separation of sandblasting material

Sand particles inside the vacuum stream can become very destructive in the long run. They can reduce the wall thickness of tubes and pipes, especially at reduction fittings and elbows, sometimes even creating holes. Sand particles should be filtered out directly where the sandblasting units are located, before they can enter the ducts of the vacuum systems. The S1 sandblasting separator is designed to accommodate up to three sandblasting units. The precipitated particles can be emptied into a simple bucket at regular intervals for disposal.





#### Planning your central vacuum system

#### CAD



Detailed planning is required to install a central vacuum system. This is the only way to provide each workspace with an adequate and uniform stream of air, regard-less of how many vacuum positi-ons are open at a time and irre-spective of their distance from the central vacuum unit. It might be helpful to call for an appointment, so our vacuum experts can visit your laboratory facilities to find out about important requirements and solutions. Detailed planning is required to

## Single suction unit

## Multi suction units

## Accessories

### **Requirements**

- You are reassured
  - about the health and environmental advantages that can only be attained with central vacuum system.
  - about the high performance and safe operation of Zubler's central vacuum system.
- You have the facilities
  - to accommodate the central unit in a separate room (e. g. a basement) or to create an appropriate installation site.
- You have reached a decision
  - whether you need separation of dust types (dirt vs. precious) metal) or whether such separation is useful given your laboratory structure.

Ask for our questionnaire if you feel these requirements are met. Complete the questionnaire and return it to us. We also need a floor plan of your existing or projected laboratory. Talk to your architect or furnisher if they can supply a computerized plan in DXF format.

### Appropriate solutions for any laboratory

It is always possible to install a central vacuum system, whether the laboratory is newly constructed or subjected to alterations, redevelopment or retrofitting. Professional and clean execution of each installation is integral to Zubler's central vacuum technology. There is nothing to hide about our equipment. However, it is both an objective and a challenge to hide it nevertheless, since our true focus is not on equipment but on offering an efficient and comfortable working environment.

#### Dental Lab Haasler Langenlois/A, 15 vac. ports

#### Laboratory redevelopment

- Vacuum pipe inside floor plate/pavement.
- Ideal routing of air streams, with the pipes taking the shortest path directly to the main connections. There is no need for paneling, as none of the pipe surfaces are visible.
- Installation is easier during laboratory construction or redesign.





#### Dental Lab Müller Dettingen, 25 vac. ports

#### Laboratory construction

- Vacuum pipes inside the floor plate (installation type offering numerous advantages).
- Configuration of workspaces as table islands.
- Requires exact planning. All installation details must be settled from the very outset of the construction works. Table groups may have to be aligned with existing floor connections.





#### Dental Lab Sobeck Lindau, 35 vac. ports

#### Laboratory redesign

- All vacuum positions integrated with furniture connected to the wall.
- "Ring circuit" concealed inside the window parapet, underneath benchtops, behind closets in the base.
- Pipe system is readily accessible and expandable.
- More complex installation, planning of extended and angulated ducts.





**Central suction system** 

#### University of Zürich Switzerland, 150 vac. ports

#### Laboratory construction

- "Double floor" approach to installation.
- Pipe system remains accessible for subsequent expansion and structural changes.
- Pipe system involving extraordinary structural requirements, since this is a school installation that needs to withstand all vacuum positions being opened at exactly the same time, which takes extra planning and dimensions.













## Accessories

To be able to use an economical and perfectly functioning suction system, it is important that all components harmonize with each other.

Therefore Zubler offers a wide range of accessories that have been created for the various applications in the dental industry.



🔺 R1200 Vacuum funnel

#### Vacuum funnels R1200 + R1250

The plug-in vacuum funnel, R1200, is used with the fixed rectangular pipe R1000. It is quiet and equipped with a large safety shield.

The vacuum funnel R1250, with a small shield can be swiveled, has

been designed for intricate work in the ceramic and precious metal department.

The gold filter can simply be exchanged if needed. It is available as an option for the vacuum funnel R1250.

> ▲ R1250 + R1300 R1200 + R1350

#### Suction hood SH-1

This suction funnel functions on every work bench without the need for permanent installation. The heavy base with a skid-proof, rubberized bottom can be used on the bench in every position. The big safety glass shield protects the technician's eyes from injuries and is particularly effective in collecting dust right where it is produced (connection  $\emptyset$  33 or 36 mm)



#### Intake system R1400

The R1400 intake unit comprises an intake funnel, an object support, and a table canal. It is mounted directly underneath the benchtop. The vacuum position will mechanically open as the suction lattice is pulled out. This system needs to be mounted

under free benchtops 20 to 40 mm in thickness. Note that its manual mechanism of vacuum position closure does not comply with certified test conditions for TÜV- or BGIAapproved vacuum systems.



#### Monomer bath R1600

By installing the monomer bath in a desk top, you obtain a comfortable stationary working place for mixing monomer and chemicals.

The mobile version R1601 can be used at different work stations by plugging it in the rectangular pipe R1000 (80x30 mm) instead of a vacuum funnel.

#### Vacuum hood R1500

Devesting and grinding creates an extreme amount of dangerous dust particles. The vacuum hood R1500, with flexible positioning arm, can be positioned according to the application.

▲ R1500

🔺 R1601

#### Vacuum plug ST02

Vacuum plugs are placed in the laboratory areas where you choose. The floor and working area are quickly and easily cleaned with cleaning accessories.



#### **Cleaning accessories**

Zubler<sup>n</sup> offers various cleaning sets as accessories. These are easily plugged into the vacuum plugs installed in the laboratory.



This allows easy cleaning of the area around the functional instruments.



**A** R1600

## AV1000

#### **Preparation** station

## For work preparation and especially dust-intensive work

Clean

Excellent dust protection via optimal placement of suction funnel, compressed air nozzle for control blowing and grinding position

Ergonomic

Large entry openings and special shaping of the arm rests allow a comfortable working position

Comfortable

Automatic vacuuming of all dust, also during controlled compressed air

#### Clearly arranged

Optimal placement of funnel, arm rests and protection glass shield provide a perfect field of view without reflection

Universal Connection to a

Connection to all vacuum systems is possible



## What is especially important for a work preparation station ?

- Complete collection of all dusts
- Collection of almost all parts which may splinter off during grinding
- Automatic vacuuming also during blowing off the dust for controlling
- Simple disposal of large pieces which may occur
- Ergonomic construction and placement of support and contact surfaces
- Noise in the collection area is reduced to a minimum



#### Cleaning

Protection shield and working block can be swiveled up out of the way. The collection tub can be taken out and easily cleaned.

#### Lighting

The work preparation station is equipped standard with a 36W Dulux energy saving bulb, giving light similar to daylight.

#### AV1000's practical and ergonomic design

A laboratory that strives for economic success in all its operations must keep a close eye on both quality and quantity. Usually this begins right in the preparatory stage for each process, where the foundation is laid for any top-quality work. This is why the design of the Zubler dust preparation box emphasizes ergonomics, cleanliness and efficiency. Stone casts are blasted clean with compressed air by actuating the integrated blast nozzle with a footswitch. You do not need to take your hands off the cast or the micromotor is an enormous time saver in preparing stone casts.

The **AV1000** grinding box comes ready for connection to external vacuum systems.





#### Technical Data AV1000

Dimension: Width x Height x Depth:	630 x 430 x 580mm
Weight:	26kg
Necessary pressure:	min. 3 bar

#### Technical Data MB1000

Dimension: Width x Height x Depth:	630 x 430 x 420mm
Weight :	18kg

## Monomer box

In the same design as the AV 1000 we also offer the Monomer station MB 1000 to capture fumes produced when working with monomer materials. It is important to use a suitable extraction system with activated charcoal filter or a central suction system that vents waste air outside the laboratory.

## P1000

#### Polishing unit

#### Polishing unit for connection to external suction systems

- Clean The basic plate is made from stainless steel and therefore resistant to aggressive acids and lye
- Quiet The integrated noise dampers reduce the noise caused by the suction air to a minimum when the polishing unit is connected to any of the Zubler<sup>®</sup> suction systems, even with very strong suction power.
- Ergonomic The working tub which gives very good protection to the working area, offers enough room to move comfortably.

## What is important for a good polishing unit ?

- Complete collection of all dusts
- Movable protection shield made of safety glass
- Optimum lighting by compact fluorescent strip bulbs
- Removable dirt tray, easy to clean and made of a soft material that will protect any falling objects.
- Intake openings offering a lownoise design, with acoustic damping of air flow and motor suspension.



P1000 for usage with zubler suction systems





#### The ergonomic concept of the **P1000** matching with the practice

The **P1000** is prepared for connection to the central suction system or a powerful external suction unit. This way, a sufficient and permanent suction efficiency is possible, even when large amounts of dust and wet pumice are involved.

It was deliberately renounced on integrating a suction system into the polishing unit. Thus, the repeated exchange of sticky filter bags is omitted, as well as the space that would have been required for the installation of a suction system. By this the complete polishing unit is smaller and can be installed even in narrow working places.

The maintenance-free short-circuited rotor offers a quiet running motor and longevity with two different speed selections.

#### **Technical Data**

Dimension: Width x Height x Depth:	630 x 460 x 490mm
Weight:	38,5kg
Power input:	550W
Speed:	1400 / 2800rpm
Noise level:	55dB(A)

## **Products for new standards**

#### Vacuum equipment

Many years of experience and continuous development have put us in a position to offer a range of vacuum solutions that will ideally meet the needs of any laboratory. Our vacuum line ranges from single-position through multiple-position units up to central systems. Small, medium and large laboratories are covered. Zubler's expertise in dental technology has proved effective around the world. Countless vacuum applications have been implemented, including complex systems with over 100 vacuum positions in schools and universities. Numerous laboratories have benefited from this expertise and are enjoying better working conditions. Our vacuum systems are complemented by low-noise intake funnels and devices for highly effective dirt absorption. All elements combined make for a perfect overall system.

#### Furnace equipment

New materials and equipment technologies are jointly developed and tested in our in-house dental laboratory. Our business philosophy is to place the same emphasis on workmanship as on design and features. VARIO PRESS® 300 is one example of the excellent outcomes of this philosophy. This ceramic furnace uses advanced sensor and control technology to ensure a uniform temperature distribution inside the firing chamber, thus offering utmost reliability of pressing and firing procedures.

#### Laboratory equipment

The P1000 polishing unit was purposefully developed without an integrated vacuum component. The unit was specifically designed for the purpose of connecting to central vacuum systems, thus optimizing vacuum performance and comfort.

#### Investment materials and pressable ceramics accessories

A variety or products are available to offer reliable results and being time savers in the hands of advanced dental technicians, such as Microstar Speed investment materials (HS-PC<sup>™</sup> for pressed ceramics, HS<sup>™</sup> for fixed crowns and partial dentures), Microstar single-use plungers for pressable ceramics, or Easy-Fix<sup>™</sup> custom firing supports.







100c

		_
Distribution:		

C E

B00092 / 03-2009



