We know how

WP Dental

**Produktkatalog
product catalogue**

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 **Special know-how, binded with amanding technologies, take care to ensure the progress in dentistry. The business we devote ourselve to, is based on recognizing the natural-science-medical connectings. Our development and combination of tailor-made problem-solvations is based hereon.**

Throughout the cooperation with universities, scientific institutions and practitioners from the area of dentistry and dental-technique, results synergical-effects which flow into our products.

According to the rules of DIN ISO 13485 all our products, courses and results get constantly proved.

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*) please ask for the catalogue of materials for the dental practise

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🇬🇧 Light curing modelling material with strong, wax-like consistency and very low contraction values. LiWa could be "hot" modelled with an electrical wax-knife (e. g. LiWaxer) or also "cold" modelled with the fingers.*

Modellations would be worked out directly onto the mastermodel, which saves not only time but also money, because there is no need for duplications. Modellations made of LiWa are polymerizing within shortest time with UV-light and provide extremely precized and accurately results.

Polymerized modellations can problem-free be taken off and prepared with rotary instruments. They dispose of excellent recovery-attitudes and high flexibility. The large assortment of preformed LiWa parts and accessories and also the different, low-priced Sets for lots of application areas leave nearly nothing to be desired.

*Assumption for modelling LiWa-material with the fingers is the wearing of adequate protection-gloves as it is advised for working with acrylates.

Indications

- For all types of metal works
- For repairs
- For crowns
- For bridges
- For implant work

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 **Isolating liquids for isolation of the master-model against the LiWa material.**

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

Indications

- Two step isolation liquids for plaster-models

Packungsart

- WP5090 LiWa Form, 12 Stk. genarbte Platten
- WP5092 LiWa Form, 30 Stk. UK-Bügel, mittel, ca. 4 mm
- WP5093 LiWa Form, 10 x 18 Stk. Molaren Klammern
- WP5094 LiWa Form, 12 Stk. Rundloch-Retention
- WP5096 LiWa Form, 30 Stk. Wachsdraht, 2 mm
- WP5096A LiWa Form, 20 Stk. Wachsdraht, 4 mm
- WP5096B LiWa Form, 30 Stk. Wachsdraht, 1 mm

Packaging

- WP5090 LiWa Form, 12 pcs. coarse-stippled plates
- WP5092 LiWa Form, 30 pcs. L-jaw bars, approx. 4 mm
- WP5093 LiWa Form, 10 x 18 pcs. Molar clasps
- WP5094 LiWa Form, 12 pcs. perforated retentions
- WP5096 LiWa Form, 30 pcs. wax-wire, 2 mm
- WP5096A LiWa Form, 20 pcs. wax-wire, 4 mm
- WP5096B LiWa Form, 30 pcs. wax-wire, 1 mm

Packungsart

- WP5095 LiWa Iso, Step I + II 2 x 10 ml

Packaging

- WP5095 LiWa Iso, Step I + II 2 x 10 ml

- WP5097 LiWa Form, 12 Stk. glatte Platten
- WP6134 LiWa Form, 5 x 20 Stk., Ringklammern
- WP6135 LiWa Form, 10 Stk. genarbte Platten, Großformat, ca. 150 x 75 mm
- WP5076 LiWa Form, 10 Stk. Retentionsgitter, runde Löcher
- WP5089 LiWa Form, 10 Stk. Retentionsgitter, quadratische Löcher

- WP5097 LiWa Form, 12 pcs. plates without structure
- WP6134 LiWa Form, 5 x 20 pcs., Ring clasps
- WP6135 LiWa Form, 10 pcs. coarse-stippled plate, large size, ca. 150 x 75 mm
- WP5076 LiWa Form, 10 pcs. Retention grid, meshes, rounded holes
- WP5089 LiWa Form, 10 pcs. Retention grid, meshes, squared holes

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UK Light curing dipping wax for producing wax caps, combined with dipping wax devices which simultaneously contain a temperature indicator.

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UK Dip Isolation based on alginate isolates plaster against plastic and is free of formaldehyde. Especially created for works with light-curing dipping-wax „LiWa Dip“

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UK Surface sealing varnish for LiWa works.

Germany LiWa Akti Bond.

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UK LiWa Akti Bond.

Indications

- Adhesive for better connectings

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- LC Tauchwachs

Indications

- LC dipping wax

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Indication

- Isolation plaster against plastic

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Indications

- For smoothing and strengthening LiWa-surfaces

Packungsart

- WP6129 LiWa Dip, transparent
2 x 40 g in Tuben
- WP5088 LiWa Dip, rot
2 x 40 g in Tuben
- WP5131 LiWa Dip, blau
2 x 40 g in Tuben

Packaging

- WP6129 LiWa Dip, transparent
2 x 40 g in tubes
- WP5088 LiWa Dip, red
2 x 40 g in tubes
- WP5131 LiWa Dip, blue
2 x 40 g in tubes

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 **LiWa Cool.**

Indications

- Cooling-spray for fast cooling of „hot“-modelled LiWa-works. Free of CFC.

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Indications

- For use as isolation against materials such as steel and plastics

Packungsart

WP5083 LiWa Universal Isolant,
10 ml

Packaging

WP5083 LiWa Universal Isolant,
10 ml

Packungsart

WP5079 LiWa Cool, Spray
1 x 200 ml

Packaging

WP5079 LiWa Cool, spray
1 x 200 ml

Packungsart

WP5084 LiWa Inhib-Clean,
250 ml

Packaging

WP5084 LiWa Inhib-Clean,
250 ml

LiWa

LiWaxer Set[®]

LiWa - Chrom - Cobalt[®]

 **LiWaxer Set.**

 **LiWaxer Set.**

 **LiWa-Chrom-Cobalt Set.**

Packungsart

WP6137 LiWa-Chrom-Cobalt Set:
LiWa I, 20 g
LiWa Iso Step I + II, 10 ml
LiWa Akti Bond, 10 ml

Packaging

WP6137 LiWa-Chrom-Cobalt Set:
LiWa I, 20 g
LiWa Iso Step I + II, 10 ml
LiWa Akti Bond, 10 ml

 **LiWa-Chrom-Cobalt Set.**

Packungsart

WP5100X LiWaxer Set:
LiWa I + II, je 20 g
LiWa Cool, 200 ml
LiWa Iso I + II, je 10 ml
Akti Bond, 10 ml
Finish Varnish, 10 ml
Universal Isolant, 10 ml
Formteil-Sortimente
LiWaxer Wachsmesser
WP5100NX LiWaxer Set
ohne Wachsmesser

Packaging

WP5100X LiWaxer Set:
LiWa I + II, 20 g each
LiWa Cool, 200 ml
LiWa Iso I + II, 10 ml each
Akti Bond, 10 ml
Finish Varnish, 10 ml
Universal Isolant, 10 ml
LiWa Form assortments
LiWaxer wax-knife
WP5100NX LiWaxer Set
without wax-knife

LiWa Form:

Molaren Klammern, 10 x 18 Stk.
eckige Lochgitter, 10 Stk.
genarbte Platten, Großformat
ca. 150 x 75 mm, 10 Stk.
UK-Bügel, ca. 4 mm, 30 Stk.
Lochretentionen, 12 Stk.
Ringklammern, 5 x 20 Stk.

LiWa Form:

Molar clasps, 10 x 18 pcs.
Retention grid meshes,
squared holes, 10 pcs.
Stippled plates, large size
approx. 150 x 75 mm, 10 pcs.
L-jaw bars, approx. 4 mm, 30 pcs.
Retentions, 12 pcs.
Ring clasps, 5 x 20 pcs.

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🇬🇧 Light curing custom tray-material, based on hybridcomposite. The highly crosslinked and always same properties allow a clean and fast work-out.

Surplus material can be reused for other works (e.g. modelling a grip for the preformed tray) without any problems.

Plaque Photo is quickly polymerized with UV- or halogenlight.

Indications

- Making of individual impression- and functional trays
- Occlusal registration
- Precized impression for part- or totalprosthetics

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Technical data

Viscosity expansion	14,0 - 19,0 mm
Daylight stability	approx. 20 min.
Flexible strength	approx. 130 MPa
Time of curing both sides (belongs to light-curing-unit)	3 - 5 min.
Plastically deformable without elastic components	
Warp resistant	

Packaging

(50 pcs. in light-protected box)

WP5001	upper-patterns, pink
WP5001A	lower-patterns, pink
WP5016	upper-patterns, blue
WP5016A	lower-patterns, blue
WP5034	upper-patterns, nature
WP5034A	lower-patterns, nature

Dentaseal[®] Pro-Dye LC[®] Pro-Dye Thinner[®]

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United Kingdom Dentaseal is a light-curing surface sealer based on multifunctional acrylates. Ready to use, for sealing and aesthetically working for dental techniques in the field of total- and partial prosthetics, individual light curing custom trays etc.

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United Kingdom Light curing one-component varnish based on multifunctional acrylate. Pro Dye LC is available in the colours: red, blue, green, yellow and nature. Throughout light-polymerization with an intensive UV/halogenlight this varnish develops a very high surface-hardness. Due to the micro-fine thickness of the layer it is able to get dry, regular and scratchresistant layer in matter of seconds.

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Indications

- Surface-sealing of individual impression trays
- as well as Part- or totalprosthesis

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Thickness of layer 10 my

Packaging

WP5515 Dentaseal, 30 ml in black brush-bottle

Packungsart

WP5518 Pro-Dye LC, blau, 30 ml in schwarzer Pinsel.
 WP5519 Pro-Dye LC, rot, 30 ml in schwarzer Pinsel.
 WP5520 Pro-Dye LC, grün, 30 ml in schwarzer Pinsel.
 WP5521 Pro-Dye LC, gelb, 30 ml in schwarzer Pinsel.
 WP5522 Pro-Dye LC, natur, 30 ml in schwarzer Pinsel.
 WP5517 Pro Dye Thinner, 10 ml Verdünnerflüssigkeit

Packaging

WP5518 Pro-Dye LC, blue, 30 ml in black brush-bottle
 WP5519 Pro-Dye LC, red, 30 ml in black brush-bottle
 WP5520 Pro-Dye LC, green, 30 ml in black brush-bottle
 WP5521 Pro-Dye LC, yellow, 30 ml in black brush-bottle
 WP5522 Pro-Dye LC, nature, 30 ml in black brush-bottle
 WP5517 Pro Dye Thinner-liquid, 10 ml

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United Kingdom *Doublident is a low-viscous floating additional curing vinyl (A)-pouring silicone for the precise pouring technique.*

Doublicates fabricated out of Doublident are tear-proof and good formable. In this case undercuts will not cause any problems.

Due to this facts Doublident guarantees a minimum of shrinkage and a long lasting storage-time, plus multiple reforming.

Doublident corresponds to the ADA Spec. No. 20.

Indications

- Pouring technique for plaster models

Technical data

Mixing time	approx. 1/2 min.
Working time (mixing time incl.)	approx. 4 min.
Curing time at 23°C	approx. 30 min.
Reforming under pressure	approx. 19 %
Shrinkage	< 0.05 %
Left over of pressure preforming	< 0.2 %
Reproduction	< 10 micron
Shore-A hardness	21/22
Colour	orange/blue

Packaging

WP5080	Doublident, 2 x 1000 ml Cat & Base, orange
WP5081	Doublident, 2 x 5000 ml Cat & Base, orange
WP5135	Doublident, blue, 2 x 1000 ml Cat & Base
WP5136	Doublident, blue, 2 x 5000 ml Cat & Base

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UK *DUR-LAB-Sil is a very hard condensation cross-linked silicone with a very high thermal stability, especially made for using in the dental-laboratory.*

Very simple handling is guaranteed because DUR-LAB-Sil does not stick, has good dimension-stabilities and is very heat-resistant.

DUR-LAB-Sil is to be used with hardener-liquid or with hardener-paste.

Indications

- Replica denture technique
- Protection of teeth when put in muffle

Technical data

Shore-A hardness + 95
Heat resistance up to 100° C

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 **Very fine spray for making contact points visible.**

Extra fine constant film thickness, fresh mint flavour and CFC-free.

Easy to be removed. Available green and red-coloured.

Indications

- For marking the contact-points at prosthetic-works

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 **Light curing, thin-flowing modelling-material for the modelation of caps and crowns.**

Due to the reinforcement with glass-particles a very high end-hardness is guaranteed after light-curing.

This high end-hardness is substantial for precized scanning for works with milling machines.

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Technical data

CFC-free
Extra fine colour pigments for exact doses

Packaging

WP5110	WP Occlusion, green, 75 ml Spray
WP5111	WP Occlusion, red, 75 ml Spray
WP5671	WP Occlusion, white, 75 ml Spray

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United Kingdom Light curing one-component-plastic for model-preparation. Useable for blocouts of undercuts, construction and completion of prepared cores and models etc. Throughout its light curing properties a fast and exact use is guaranteed.

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
New Pontics made of plastics with benefits like:

- Delivers accurately fitting castings and scan objects through the form-stability
- Residue-free and smokeless burnings
- Usable for gypsum-bonded investment
- Moldable (even in the length) after heating
- Possible binding with modeling-materials
- Clean working, because the product is non greasy and non-recovery
- Easy milling – therefore easy connecting especially for implant-works
- Biostatic binding-design
- Ideal for copy-milling

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Packungsart

WP4060 Extraclean, 1 l
WP4061 Extraclean, 5 l

 **Concentrate for easy cleaning of impression trays and instruments. Used in a 5 % solution it removes residues from silicones, alginates, cements easily from trays and dental-instruments.**

Indications

- Cleaning

Technical data

Concentration 5 % / liter
UBA no.: 23660001
(UBA = environmental office of Federal Republic of Germany)

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 **Wapol is a powerful, non-emulgating special wax-remover. Wapol guarantees a thorough, residue-free and intensive elimination/cleaning of any kind of hot-wax devices or wax-melting-pots. Its concentrated form allows a very economical use.**

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 **Calsol is a pH-neutral, plaster and stone-remover.**

If used according to the instructions it removes easily the dental plaster from prosthesis and impression trays.

Calsol is ready for use. If used in ultrasonic it dilutes even larger parts of dental plaster within shortest time.

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WP5588 Wapol, 1.000 ml
WP5589 Wapol, 5.000 ml

Indications

- Cleaning of hot-wax-baths and wax-melting-pots

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WP5590 Calsol, 1.000 ml
WP5591 Calsol, 5.000 ml

Indications

- Cleaning of dental plaster contaminated prosthesis and equipment

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
Packungsart

WP6018 LiWa Light, 1 Stck.

UV-light curing unit, especially designed for LiWa-works. The built-in ventilator cools the LiWa-material while polymerizing it. Therefore LiWa does not run off the model and a perfect polymerizing-result is guaranteed. LiWa Light could be timed with freely choosable times (secondly).

Technical data

Operating voltage	230 V/115V
Power frequency	AC 50/60 Hz
Fuse link	1A
Spectral range	UVA
Light source	2 PL-S 9W

Weight 2,8 kg
Size over all 8,8x13.5x28 cm
Self-starting ventilator for cooling of LiWa[®].

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Packungsart

WP6002 Highlight Easy
+ 4 x 9 W UV-Lampen
WP6006 Ersatzlampe

UVA-light curing unit (tunnel-version) for polymerization of light curing custom tray material. This unit contains a switch with two different preselected curing-times (120 sec. / 180 sec.).

Technical data

36 Watt =	4 x 9 Watt UVA-bulbs
Times =	120 sec. / 180 sec.

Packungsart

WP6002 Highlight Easy
+ 4 x 9 W UVA-bulbs
WP6006 replacement bulb

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 **Compact electrical waxknife for LiWa works with infinitely variable temperature-adjustment. It contains an ergonomically formed handpiece with easily to change modelling-tips and a slightly preheating-phase of only 20-30 seconds with 7.5 W.**

Technical data

Measures (L/W/H) 100 x 65 x 55 mm
Voltage 220-240 V/50-60 Hz
Plunger/Handpiece 12 V / 7.5 W
Variable temperature-adjustment 50° C - 350° C
Heating-control LED
Optional accessories:
3 modeling-tips (blade, spatula, knife)

Packaging

WP6110 1 LiWaxer modeling-station
+ modeling-awl
+ plunger-depot
WP5105 Replacement-tip (Awl)
WP5106 Replacement-tip (Knife)
WP5107 Replacement-tip (Blade)
WP5108 Replacement-tip (Spatula)
WP5109 Replacement Soldering Iron

We know how

Willmann & Pein GmbH

Schusterring 35
D-25355 Barmstedt/Hamburg
Phone: +49-4123-9228-0
Fax: +49-4123-9228-49
info@wp-dental.de
www.wp-dental.de