


EVE TWIST



DATA SHEET

EVE TWIST





Item No.: KT9601 | Edition: 04/2013




EVE TWIST | Metal processing

EVE TWIST reduces the need for the dental technician to resort to the traditional high shine polishing with brushes and paste. EVE TWIST is based on the concept of thin discs. Depending on the type of application the dental technician needs he can mount two or more layers on top of the other.


Due to the innovative radial shape and the mounting options the processed surface is not altered during the final high-shine polishing process.



Combination of flat surfaces = Geometry Lens



Combination of arched surfaces = Geometry Wheel



extra coarse
 coarse
 medium
 fine

EVE TWIST

Pack sizes: 10, 100 pieces

rec. RPM:

7.000 – 12.000 min

Item::
Item No.:
Sizes:
Grit:

ET-514	ET-614	ET-714	ET-814
4514	4614	4714	4814
14 x 1,2	14 x 1,2	14 x 1,2	14 x 1,2
extra coarse	coarse	medium	fine

Pack sizes: 10, 100 pieces

rec. RPM:

7.000 – 12.000 min

Item::
Item No.:
Sizes:
Grit:

ET-517	ET-617	ET-717	ET-817
4517	4617	4717	4817
17 x 1,2	17 x 1,2	17 x 1,2	17 x 1,2
extra coarse	coarse	medium	fine

Pack sizes: 10, 100 pieces

rec. RPM:

7.000 – 12.000 min

Item::
Item No.:
Sizes:
Grit:

ET-526	ET-626	ET-726	ET-826
4526	4626	4726	4826
26 x 1,2	26 x 1,2	26 x 1,2	26 x 1,2
extra coarse	coarse	medium	fine