



W o r l d L e a d e r i n S p h e r i c a l A t t a c h m e n t s


PRICE LIST

for Dentist and Dental Laboratories



[2012]

w w w . r h e i n 8 3 . c o m





CASTABLE



Ref. 092ECQ

OT EQUATOR CASTABLE

KIT CONTAINS:

- 2 SINGLE MALES
- 2 STAINLESS STEEL HOUSINGS
- 4 RETENTIVE CAPS
(2 soft, 2 standard)

€ 50.00

IMPLANTS



Ref. 130

OT EQUATOR for IMPLANTS

KIT CONTAINS:

- 1 TITANIUM ABUTMENT
- 1 STAINLESS STEEL HOUSINGS
- 1 PROTECTIVE DISK
- 4 RETENTIVE CAPS
(1 extra-soft, 1 soft, 1 standard, 1 strong)

€ 115.00



Ref. 030

- 1 TITANIUM ABUTMENT

€ 80.00

BARS



Ref. 160EQB

OT EQUATOR threaded sleeve system for bonding

KIT CONTAINS:

- 2 EQUATOR THREADED for titanium sleeve (1.6 mm thread)
- 2 THREADED SLEEVES (1.6 mm thread)
- 2 STAINLESS STEEL HOUSINGS
- 2 EQUATOR SPACERS
- 8 RETENTIVE CAPS
(2 extra-soft, 2 soft, 2 standard, 2 processing)

€ 125.00



Ref. 039SFE2

- 1 THREADED EQUATOR (direct into 2 mm thread)

€ 25.00

Passive bar connections system "ELASTIC SEEGER"



Ref. 158EQS

ELASTIC SEEGER

KIT CONTAINS:

- 2 CASTABLE ABUTMENTS FOR SEEGER
- 3 RED PLASTIC SEEGER FOR LABORATORY
- 3 RETENTIVE ELASTIC FOR LOCKING THE ABUTMENTS
- 2 TOP LOCKING TITANIUM SCREWS

€ 45.00

SPARE PARTS:



Ref. 192ECE

CAPS ASSORTMENT KIT

KIT CONTAINS:

- 1 STAINLESS STEEL HOUSINGS
- 1 PROCESSING CAP
- 4 RETENTIVE CAPS
(1 extra-soft, 1 soft, 1 standard, 1 strong)

€ 38.00



Ref. 151SS

• 2 CASTABLE MALES

€ 20.00



Ref. 140CEV

• 4 RETENTIVE CAPS - violet "STRONG" (2.7kg)

€ 28.00



Ref. 140CET

• 4 RETENTIVE CAPS - white "STANDARD" (1.8kg)

€ 28.00



Ref. 140CER

• 4 RETENTIVE CAPS - pink "SOFT" (1.2kg)

€ 28.00



Ref. 140CEG

• 4 RETENTIVE CAPS - yellow "EXTRA-SOFT" (0.6kg)

€ 28.00



Ref. 140CEN

• 4 PROCESSING CAPS - black for lab

€ 14.00



Ref. 141CAE

• 2 STAINLESS STEEL HOUSINGS

€ 20.00



Ref. 339SFE

• 1 EQUATOR THREADED 1.6mm (for titanium sleeve)

€ 25.00



Ref. 239GSF

• 1 TITANIUM SLEEVE (1.6mm thread)

€ 24.00



Ref. 339SSE

• 1 EQUATOR SPACER

€ 10.00



Ref. 158TT

• 1 TITANIUM LOCKING SCREW

€ 14.00

ACCESSORIES:



Ref. 144MTE

• 2 IMPRESSION COPING
for individual tray

€ 36.00



Ref. 144AE

• 2 LABORATORY ANALOG

€ 38.00



Ref. 044CAIN

• 2 IMPRESSION COPINGS NORMAL SIZE - OT EQUATOR
with interchangeable cap for pick up impression

€ 25.00

TOOLS:



Ref. 74AC01

• 1 PARALLELOMETER KEY "NORMAL"

€ 22.00



Ref. 774CQ

• 1 OT-EQUATOR SQUARE SCREW DRIVER
for implant abutment (square 1.25mm)

€ 39.00



Ref. 760CE

• 1 SQUARE DRIVER CONNECTOR for torque
controller contra angle (square 1.25mm)

€ 39.00



Ref. 124ICP

• 1 BLUE PLASTIC "MULTIUSE" INSERTION TOOL
(used together with metal insertors)

€ 15.00



Ref. 185IAC

• 1 METAL INSERTION TOOL FOR CAPS
(use together with 124ICP BLU HANDLE)

€ 10.00



Ref. 185SIS

• 1 METAL INSERTION TOOL FOR SEEGER
(use together with 124ICP BLU HANDLE)

€ 12.00



Ref. 191ECS

• 1 METAL EXTRACTOR TOOL FOR CAPS
(use together with 124ICP BLU HANDLE)

€ 16.00

NEW

OT UNILATERAL

EXTRA CORONAL ATTACHMENTS WITH COMBINED FUNCTION



Ref. 098UNL

OT UNILATERAL + CONNECTOR

KIT CONTAINS:

- 1 UNILATERAL CASTABLE ATTACHMENT
- 1 UNI-BOX CASTABLE
- 1 PLASTIC POSITIONER MICRO SIZE
- 1 PINK CAP MICRO SIZE
- 1 STRATEGY YELLOW CAP (for duplication technique)
- 1 BLACK CAP MICRO SIZE (for laboratory use only)
- 1 STRATEGY BLACK CAP (for laboratory use only)
- 1 CONNECTOR

€ 35.00



Ref. 156BCU

- 1 UNILATERAL CASTABLE ATTACHMENT

€ 12.00



Ref. 156CCU

- 1 UNI-BOX CASTABLE

€ 12.00

NEW

OT CAP IMPRESSION COPINGS

NORMAL - MICRO



Ref. 044CAIN

- 2 IMPRESSION COPINGS NORMAL SIZE
with interchangeable cap

€ 25.00



Ref. 044CAIM

- 2 IMPRESSION COPINGS MICRO SIZE
with interchangeable cap

€ 25.00

SINGLE THREADED SPHERES

with THREADED SLEEVE TO BE GLUED NORMAL - MICRO



Ref. 139KSFN

KIT SINGLE THREADED SPHERE NORMAL SIZE + TiN with sleeve to be BONDED

KIT CONTAINS:

- 2 TITANIUM NORMAL SINGLE THREADED SPHERE
(HEX mm. 1.3 THREAD mm. 1.6)
- 2 TITANIUM BOND-IN THREADED SLEEVE
- 2 S/STEEL SPACER for the THREADED SLEEVE (Normal sphere)

€ 88.00



Ref. 139KSFM

KIT SINGLE THREADED SPHERE MICRO SIZE + TiN with sleeve to be BONDED

KIT CONTAINS:

- 2 TITANIUM MICRO SINGLE THREADED SPHERE
(HEX mm. 0.9 THREAD mm. 1.6)
- 2 TITANIUM BOND-IN THREADED SLEEVE
- 2 S/STEEL SPACER for the THREADED SLEEVE (Micro sphere)

€ 88.00

SPARE PARTS:



Ref. 239SFN

- 1 SINGLE THREADED SPHERE NORMAL SIZE + TiN
for threaded sleeve
HEX mm. 1.3 THREAD mm. 1.6

€ 20.00



Ref. 239SFM

- 1 SINGLE THREADED SPHERE MICRO SIZE + TiN
for threaded sleeve
HEX mm. 0.9 THREAD mm. 1.6

€ 20.00



Ref. 239GSF

- 1 TITANIUM BOND-IN THREADED SLEEVE
for Normal and Micro - mm. 1.6

€ 24.00



Ref. 239SSN

- 1 STAINLESS STEEL SPACER for THREADED SLEEVE
(Normal sphere)

€ 10.00



Ref. 239SSM

- 1 STAINLESS STEEL SPACER for THREADED SLEEVE
(Micro sphere)

€ 10.00

NEW

SINGLE THREADED SPHERES

For Cad/Cam or Self Tapping 2.0mm Thread - Normal and Micro Size



Ref. 039SFN2

- **1 SINGLE THREADED SPHERE NORMAL SIZE + TiN**
HEX 1.3mm, THREAD 2.0mm

€ 20.00



Ref. 039SFM2

- **1 SINGLE THREADED SPHERE MICRO SIZE + TiN**
HEX 0.9mm, THREAD 2.0mm

€ 20.00



Ref. 760CBR

- **HEX DRIVER FOR THREADED SPHERE NORMAL**
HEX 1.3mm

€ 39.00



Ref. 772CSF

- **HEX DRIVER FOR THREADED SPHERE MICRO**
HEX 0.9mm

€ 39.00

OT CAP & OT BOX - OT STRATEGY - OT BAR - OT VERTICAL ATTACHMENTS OT TECNO - TITANIUM SINGLE SPHERE

MASTER STARTER KIT FOR LABORATORY:



Ref. 005SKLUS

KIT CONTAINS:

OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (normal and micro size)
- 4 castable spheres (normal and micro size)
- 2 titanium single spheres (normal and micro size) + sphere holder + 2 sliding base
- 4 castable OT CAP bars (normal and micro size)
- 2 complete castable OT CAP TECNO bars (normal and micro size)
- 4 castable OT BOX CLASSIC bars (upper + lower part) normal and micro size
- 2 castable OT BOX SPECIAL bars (normal and micro size)
- 6 castable OT BOX connectors
- 4 castable OT BOX MONO housings (normal and micro size)
- 4 plastic positioner rings (normal and micro size)
- 4 OT CAP retentive caps (normal and micro size)
- 16 WHITE- standard retention - PINK soft retention - YELLOW extra soft retention
GREEN very elastic retention-
- 4 TITAN CAP titanium caps (normal and micro size)
- 4 BLACK CAPS- only for lab use - (normal and micro size)
- 8 S/STEEL housings
- 4 ANALOGUES for plaster models (normal and micro size)
- 4 protective disks
- 4 aluminium spacers (normal and micro size)

OT STRATEGY:

- 4 castable spheres (2 standard and 2 high)
- 2 castable STEADY
- 6 retentive caps: 2 white - 2 pink - 2 yellow for housing
- 2 black caps for housing only for lab use
- 2 s/steel housings
- 2 plastic positioner rings
- 6 retentive caps for duplication technique:
WHITE - standard retention - PINK soft retention - YELLOW extra soft retention
- 2 black caps for duplication technique only for lab use
- 2 ANALOGUES for plaster models

OT BAR MULTIUSE:

- 2 castable bar
- 1 gingival connector
- 2 positioner clip A
- 2 positioner clip B
- 4 castable housings
- 4 clips - 2 PINK - soft retention - 2 YELLOW medium retention

OT VERTICAL:

- 2 castable male
- 2 castable STEADY
- 4 clips - 2 WHITE- standard retention - 2 GREEN soft retention
- 2 castable parallelometer keys
- 2 ceramic bars

TOOLS:

- 1 Pincers
- 2 Reamer tools for caps (normal and micro size)
- 1 Reamer tools for OT STRATEGY caps
- 2 Parallelometer keys for OT CAP (normal and micro size)
- 1 Parallelometer key for OT STRATEGY size
- 2 Parallelometer keys for OT CAP TECNO (normal and micro size)
- 1 Parallelometer key for OT BAR MULTIUSE
- 1 Blue Plastic "multiuse" insertion tool
(for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- 1 OT STRATEGY insertion tool
- 1 OT BAR clips insertion tool
- 1 OT VERTICAL clips insertion tool

€ 400.00

OT CAP & OT BOX - OT STRATEGY - OT BAR - OT VERTICAL ATTACHMENTS

BASIC LABORATORY KIT:



Ref. 005SKLBUS

NEW

KIT CONTAINS:

OT CAP - OT BOX SYSTEM:

- 16 castable assorted pivots (normal and micro size)
- 2 castable spheres normal
- 2 castable spheres micro
- 2 castable OT CAP bars (normal and micro size)
- 2 castable OT BOX CLASSIC bars (upper + lower part) normal and micro size
- 1 castable OT BOX SPECIAL bars (normal and micro size)
- 6 castable OT BOX connectors
- 4 castable OT BOX MONO housings (normal and micro size)
- 8 plastic positioner rings (normal and micro size)
- 28 OT CAP assorted retentive caps (normal and micro size)
WHITE- standard retention - PINK soft retention - YELLOW extra soft retention
GREEN very elastic retention
- 9 BLACK CAPS- only for lab use - (normal and micro size)
- 4 S/STEEL housings (normal and micro size)

OT STRATEGY:

- 4 castable spheres (2 standard and 2 high)
- 2 castable STEADY
- 8 retentive caps: 2 white - 2 pink - 2 yellow - 2 black for housing
- 2 s/steel housings
- 2 plastic positioner rings
- 6 retentive caps for duplication technique:
WHITE - standard retention - PINK soft retention - YELLOW extra soft retention
- 2 black caps for duplication technique only for lab use

OT EQUATOR:

- 2 single males
- 2 stainless steel housings
- 4 retentive caps (2 soft, 2 standard)
- 2 processing black caps

OT UNILATERAL:

- 1 unilateral castable attachment
- 1 uni-box castable
- 1 plastic positioner micro size
- 2 ot cap micro size: 1 pink - 1 black processing cap
- 2 ot strategy (for duplication technique): 1 pink - 1 black processing cap
- 1 connector

OT BAR MULTIUSE:

- 1 castable bar
- 1 gengival connector
- 4 positioner clip (A - B)
- 2 castable housings
- 4 clips - 2 PINK - soft retention - 2 YELLOW medium retention

OT VERTICAL:

- 2 castable male
- 2 castable STEADY
- 4 clips - 2 WHITE- standard retention - 2 GREEN soft retention
- 2 castable parallelometer keys
- 2 ceramic bars

TOOLS:

- 1 Pincers
- 1 Parallelometer keys for OT CAP normal
- 1 Parallelometer keys for OT CAP micro
- 1 Parallelometer key for OT STRATEGY size
- 1 Parallelometer key for OT BAR MULTIUSE
- 1 Blue Plastic "multiuse" insertion tool
(for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- 1 OT STRATEGY insertion tool
- 1 OT BAR clips insertion tool
- 1 OT VERTICAL clips insertion tool

€ 150.00

ATTACHMENTS OT CAP SYSTEM



Ref. 038STF
NORMAL SIZE

TITANIUM FLEX SINGLE SPHERE NORMAL SIZE + TiN
FOR LASER SOLDERING - FOR BONDING

KIT CONTAINS:

- 2 SINGLE TITANIUM FLEX SPHERES
- 2 PINK CAPS - SOFT RETENTION
- 1 SPHERE HOLDER
- 2 CASTABLE SLIDING BASES

€ 75.00



Ref. 038STN
NORMAL SIZE

TITANIUM SINGLE SPHERE NORMAL SIZE + TiN
FOR LASER SOLDERING - FOR BONDING

KIT CONTAINS:

- 2 SINGLE TITANIUM SPHERES
- 2 PINK CAPS - SOFT RETENTION
- 1 SPHERE HOLDER
- 2 CASTABLE SLIDING BASES

€ 50.00



Ref. 038STM
MICRO SIZE

TITANIUM SINGLE SPHERE MICRO SIZE + TiN
FOR LASER SOLDERING - FOR BONDING

KIT CONTAINS:

- 2 SINGLE TITANIUM SPHERES
- 2 PINK CAPS - SOFT RETENTION
- 1 SPHERE HOLDER
- 2 CASTABLE SLIDING BASES

€ 50.00



Ref. 092CAN

OT CAP NORMAL SIZE

KIT CONTAINS:

- 4 SINGLE SPHERES
- 4 PINK RETENTIVE CAPS
- 4 STAINLESS STEEL HOUSINGS
- 4 PLASTIC POSITIONERS

€ 49.00



Ref. 092CAM

OT CAP MICRO SIZE

KIT CONTAINS:

- 4 SINGLE SPHERES
- 4 PINK RETENTIVE CAPS
- 4 STAINLESS STEEL HOUSINGS
- 4 PLASTIC POSITIONERS

€ 49.00



Ref. 196BCN

“ECONOMIC” OT CAP NORMAL SIZE

KIT CONTAINS:

- 1 CASTABLE PRESECTIONED BAR
- 1 CASTABLE PRESECTIONED BEVELED BAR
- 4 WHITE RETENTIVE CAPS (Standard retention)

€ 39.00



Ref 197BCM

“ECONOMIC” OT CAP MICRO SIZE

KIT CONTAINS:

- 1 CASTABLE PRESECTIONED BAR
- 1 CASTABLE PRESECTIONED BEVELED BAR
- 4 WHITE RETENTIVE CAPS (Standard retention)

€ 39.00



Ref. 099BSN

OT CAP & Mono OT BOX FOR FRAME NORMAL SIZE

KIT CONTAINS:

- 2 CASTABLE PRESECTIONED BARS (1 straight + 1 beveled)
- 4 WHITE RETENTIVE CAPS
- 4 CASTABLE MONO OT BOX
- 4 PLASTIC POSITIONERS

€ 64.00



Ref. 099BSM

OT CAP & Mono OT BOX FOR FRAME MICRO SIZE

KIT CONTAINS:

- 2 CASTABLE PRESECTIONED BARS (1 straight + 1 beveled)
- 4 WHITE RETENTIVE CAPS
- 4 CASTABLE MONO OT BOX
- 4 PLASTIC POSITIONERS

€ 64.00



Ref. 090TCN

OT CAP TECNO - TITAN CAP NORMAL SIZE

KIT CONTAINS:

- 1 OT CAP TECNO PRESECTIONED BAR
- 2 CONCAVE TITANIUM SPHERES - TIN COATED
- 2 TITAN CAPS
- 1 INSERTING TOOL

€ 110.00



Ref. 090TCM

OT CAP TECNO - TITAN CAP MICRO SIZE

KIT CONTAINS:

- 1 OT CAP TECNO PRESECTIONED BAR
- 2 CONCAVE TITANIUM SPHERES - TIN COATED
- 2 TITAN CAPS
- 1 INSERTING TOOL

€ 110.00



Ref. 058BSN

OT BOX SPECIAL NORMAL SIZE + CONNECTORS

KIT CONTAINS:

- 2 OT BOX SPECIAL BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

€ 37.00



Ref. 058BSM

OT BOX SPECIAL MICRO SIZE + CONNECTORS

KIT CONTAINS:

- 2 OT BOX SPECIAL BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

€ 37.00



Ref. 153BCN

OT BOX CLASSIC NORMAL SIZE + CONNECTORS

KIT CONTAINS:

- 2 UPPER BARS
- 2 LOWER BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

€ 34.00



Ref. 153BCM

OT BOX CLASSIC MICRO SIZE + CONNECTORS

KIT CONTAINS:

- 2 UPPER BARS
- 2 LOWER BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

€ 34.00

OT BOX LARGE

OVERSIZED CASTABLE HOUSING FOR POSITIONING CAPS CHAIRSIDE



Ref. 253BLN

OT BOX LARGE NORMAL & CONNECTORS

KIT CONTAINS:

- 2 OT BOX LARGE BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

€ 34.00



Ref. 253BLM

OT BOX LARGE MICRO & CONNECTORS

KIT CONTAINS:

- 2 OT BOX LARGE BARS
- 4 PLASTIC POSITIONERS
- 4 CONNECTORS

€ 34.00

OT STRATEGY ATTACHMENTS



Ref. 098SSS

OT STRATEGY

CAPS FOR DUPLICATION TECHNIQUE

KIT CONTAINS:

- 4 CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)
- 2 CASTABLE STEADY
- 4 RETENTIVE CAPS FOR DUPLICATION TECHNIQUE

€ 49.00



Ref. 098CAL

OT STRATEGY

CAPS FOR STAINLESS STEEL HOUSING TECHNIQUE

KIT CONTAINS:

- 4 CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)
- 2 CASTABLE STEADY
- 4 STAINLESS STEEL HOUSINGS
- 2 PLASTIC POSITIONERS
- 4 RETENTIVE CAPS FOR HOUSING TECHNIQUE

€ 59.00

OT BAR MULTIUSE ATTACHMENTS



Ref. 021OBM

OT BAR MULTIUSE & CONNECTORS

KIT CONTAINS:

- 2 CASTABLE BARS
- 4 POSITIONER CLIP A
- 4 POSITIONER CLIP B
- 4 CASTABLE BOXES
- 4 RETENTIVE PINK CLIPS
- 4 RETENTIVE YELLOW CLIPS
- 2 CONNECTORS
- 1 GINGIVAL CONNECTOR

€ 92.00

OT VERTICAL ATTACHMENTS



Ref. 071OBV

OT VERTICAL


KIT CONTAINS:

- 4 CASTABLE MALES
- 2 CASTABLE STEADY
- 4 RETENTIVE WHITE CLIPS
- 2 RETENTIVE GREEN CLIPS
- 4 CERAMIC BARS
- 4 CASTABLE PARALLELOMETER KEYS + PIN

€ 50.00

OT REVERSE 3

EXTRA-CORONAL ATTACHMENTS: castable bars and retentive females in titanium + TiN to be bonded



Ref. 095PCR

OT REVERSE 3

KIT CONTAINS:

- 2 Castable bars
- 2 Retentive females in Titanium+TiN to be bonded
- 2 Retentive males attachment Titanium + nylon
- 2 Positioning rings for OT REVERSE 3
- 2 OT box mono

€ 120.00

OVERDENTURE ATTACHMENTS



Ref. 011SPL

S.P.L. PIVOTS - COPING COVER PROMOTIONAL KIT FOR CLINIC

KIT CONTAINS:

- 24 CASTABLE ASSORTED PIVOTS MICRO SIZE
- 40 CASTABLE ASSORTED PIVOTS NORMAL SIZE
- 12 TITANIUM MICRO PIVOTS (adapted for receiving COPING COVER)
- 12 TITANIUM NORMAL PIVOTS (adapted for receiving COPING COVER)
- 16 NORMAL CAPS (WHITE - standard retention - PINK soft retention)
- 16 MICRO CAPS (WHITE - standard retention - PINK soft retention)
- 20 PROTECTIVE DISKS
- 2 REGULATING TOOLS (1 normal and 1 micro size)
- 1 BLUE PLASTIC "MULTIUSE" INSERTION TOOL (for steel tool and for inserting the caps in the mouth)
- 1 OT CAP insertion tool
- 4 NORMAL STAINLESS STEEL HOUSING
- 4 MICRO STAINLESS STEEL HOUSING
- 12 DIRECTIONAL RINGS (grey 0° - green 7° - red 14°)
- 1 PINCERS

€ 310.00



Ref. 033PSF

PIVOT FLEX TiN COATED

TITANIUM PIVOT WITH SWIVEL SPHERE - Normal size
Adapted for COPING COVER

KIT CONTAINS:

- 1 TITANIUM PIVOT WITH SWIVEL SPHERE (adapted for COPING COVER)
- 1 STAINLESS STEEL HOUSING
- 2 PINK CAPS Normal Size - Soft retention
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

€ 48.00

Available for lengths
Normal sphere ø 2,5
Micro sphere ø 1,8



Ref. 036PTN

TITANIUM PIVOTS KIT NORMAL SIZE

Adapted for COPING COVER

KIT CONTAINS:

- 5 TITANIUM PIVOTS (Adapted for COPING COVER)

€ 39.00

Ref. 036PTM

TITANIUM PIVOTS KIT MICRO SIZE

Adapted for COPING COVER

KIT CONTAINS:

- 5 TITANIUM PIVOTS (Adapted for COPING COVER)

€ 39.00

OT REVERSE 3

OVERDENTURE POSTS: for roots, available in one standard size

STANDARD LENGTH



Ref. 034PRK

KIT CONTAINS:

- 2 Root posts in Titanium + TiN
- 2 Retentive male attachment Titanium + nylon
- 2 Plastic manual tools
- 2 Stainless steel housings
- 2 Pink protective disks

€ 120.00



Ref. 037CCL

COPING COVER FOR OVERCASTING

KIT CONTAINS:

- 1 ASSEMBLING RING FOR PIVOTS
- 2 SPRINGS

€ 35.00

Ref. 010PSP

CASTABLE PIVOTS FOR CANAL IMPRESSIONS

KIT CONTAINS:

- 20 CASTABLE PIVOTS NORMAL SIZE

€ 24.00

Ref. 012PSM

CASTABLE PIVOTS FOR CANAL IMPRESSIONS

KIT CONTAINS:

- 20 CASTABLE PIVOTS MICRO SIZE

€ 26.00



Ref. A01MOG

MOOSER BUR

FOR PIVOTS OF 7 - 9 - 10mm.

€ 15.00



Ref. A03MOB

MOOSER BUR

FOR PIVOTS OF 12 - 14mm.

€ 15.00



Ref. 121FD

OT REVERSE 3 DIAMOND SHAPING BUR FOR ROOTS

€ 25.00

IMPLANTOLOGY TITANIUM ATTACHMENTS FOR OVERDENTURE
SPHERO FLEX - BLOCK SYSTEM:



SPHERO FLEX WITH TiN COATED

Ref. 109

KIT CONTAINS:

- 1 TITANIUM ABUTMENT WITH SWIVEL SPHERE (SPHERE Ø 2.5 mm.)
- 2 PINK CAPS (Soft retention)
- 1 STAINLESS STEEL HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

€ 103.00



SPHERO BLOCK NORMO WITH TiN COATED

Ref. 002

KIT CONTAINS:

- 1 TITANIUM ABUTMENT (SPHERE Ø mm. 2.5)
- 2 PINK CAPS (Soft retention)
- 1 STAINLESS STEEL HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

€ 85.00



SPHERO BLOCK MICRO WITH TiN COATED

Ref. 003

KIT CONTAINS:

- 1 TITANIUM ABUTMENT (SPHERE Ø mm. 1.8)
- 2 PINK CAPS (Soft retention)
- 1 STAINLESS STEEL HOUSING
- 1 PROTECTIVE DISK
- 3 DIRECTIONAL RINGS

€ 85.00




OT REVERSE 3 WITH TiN COATED

Ref. 120

- 1 OT REVERSE 3 (titanium+TiN) implant abutment
- 1 RETENTIVE MALE ATTACHMENT (titanium + nylon)
- 1 PLASTIC MANUAL TOOL
- 1 STAINLESS STEEL HOUSING
- 1 PINK PROTECTIVE DISK

€ 90.00

The SPHERO FLEX, SPHERO BLOCK and OT REVERSE attachments are available for the major implant brands.
Upon request we can customize them for any kind of implants.
When ordering, please always indicate: implant brand, implant diameter and cuff height.
FOR SPHERO FLEX ABUTMENT CUFF AVAILABLE FROM 1 TO 7 MM.



AUXILIARY TOOLS: SPHERO FLEX - BLOCK; OT REVERSE 3



Ref. 771CEF

KEY FOR SPHERO FLEX AND SPHERO BLOCK (Normal - Micro size)
HEX mm. 2.3

€ 48.00



Ref. 760CBM

HEX SCREWDRIVER FOR TORQUE CONTROLLER
CONTRA ANGLE

€ 39.00



Ref. 773CAR

OT REVERSE 3 ABUTMENT WRENCH

€ 60.00

ANTI-UNSCREWING SYSTEM WITH AN ELASTIC DOWEL APPLICABLE TO

SPHERO FLEX-BLOCK - OT REVERSE AND DIFFERENT SCREWS. THE ELASTIC DOWEL IS APPLICABLE ONLY UPON REQUEST TO ALL SCREWS WITH A NOT INFERIOR M. 1.8 to screws - Sphero Flex - Sphero Block - OT Reverse attachments, etc.

Ref. TE

Extra charge for applying the dowel



€ 20.00

AUXILIARY TOOLS SPHERO FLEX - BLOCK



Ref. 100AD

DIRECTIONAL RINGS FOR SPHERO FLEX AND BLOCK
KIT CONTAINS:
• 3 COLOURED RINGS WITH INCLINATION OF 0° - 7° - 14°

€ 10.00



Ref. 044ABN

SPHERO FLEX - BLOCK NORMO ANALOG KIT:
KIT CONTAINS:
• 1 SPHERO BLOCK-FLEX NORMO ANALOG
• 3 DIRECTIONAL RINGS (0° - 7° - 14°)
• 1 BLACK NORMO CAP FOR LABORATORY USE

€ 18.00



Ref. 044ABM

SPHERO BLOCK MICRO ANALOG KIT:
KIT CONTAINS:
• 1 SPHERO BLOCK-MICRO ANALOG
• 3 DIRECTIONAL RINGS (0° - 7° - 14°)
• 1 BLACK MICRO CAP FOR LABORATORY USE

€ 18.00

CUFF HEIGHT MEASURER


This tool is used to determine the cuff height of implants with internal or external hexagon.

NEW



Ref. 008MBG

KIT CONTAINS:
• 1 CUFF HEIGHT SLIDER GAUGE
• 1 CUFF HEIGHT FIXED ROD GAUGE
• 1 SILICON RING DISPENSER
• 20 SILICON RINGS



€ 80.00

OT LOCK LOCKING PIN TITANIUM AND CASTABLE



OT LOCK KIT CASTABLE

Ref. 880CLT

KIT CONTAINS:
• 1 LOCKING PIN SET
• 1 HOUSING SHAPER
• 1 CERAMIC PIN

€ 70.00



OT LOCK KIT CASTABLE

Ref. 880CLC

KIT CONTAINS:
• 1 LOCKING PIN SET
• 2 CERAMIC PIN
• 1 CERAMIC SPACER AND CASTABLE RING
• 1 HOUSING SHAPER

€ 65.00



OT LOCK KIT TITANIUM ADJUSTABLE

Ref. 880CLR

KIT CONTAINS:
• 1 LOCKING PIN SET
• 1 HOUSING SHAPER
• 1 CERAMIC PIN
• 9 CASTABLE SPACER RINGS

€ 70.00

OT LOCK ACCESSORIES



Ref. 882CG

• 1 CONICAL GUIDE

€ 10.00



Ref. 882CAS

• 1 UNLOCKING TOOL

€ 15.00

SPARE PARTS: OT CAP SYSTEM








NORMAL Ø 2.5 mm:

	Ref. 055SCN	• 4 CASTABLE SINGLE SPHERES NORMAL SIZE	€ 19.00
	Ref. 151BAN	• 2 CASTABLE PRESECTIONED BEVELED BARS NORMAL SIZE	€ 19.00
	Ref. 150BPN	• 2 CASTABLE PRESECTIONED BARS NORMAL SIZE	€ 19.00
	Ref. 048CAN	• 6 SILVER CAPS NORMAL SIZE (EXTRA-RESILIENT) ELASTIC AND GUMMY	€ 26.00
	Ref. 048CON	• 6 GOLD CAPS NORMAL SIZE (EXTRA-RESILIENT) SLIGHTLY ELASTIC	€ 26.00
	Ref. 049PCN	• 6 GREEN CAPS NORMAL SIZE (VERY ELASTIC RETENTION)	€ 26.00
	Ref. 060CRN AY	• 6 YELLOW CAPS NORMAL SIZE (EXTRA SOFT RETENTION)	€ 26.00
	Ref. 040CRN SN	• 6 PINK CAPS NORMAL SIZE (SOFT RETENTION)	€ 26.00
	Ref. 040CRN	• 6 WHITE CAPS NORMAL SIZE (STANDARD RETENTION)	€ 26.00
	Ref. 040TCN	• 2 TITAN CAPS NORMAL SIZE	€ 36.00
	Ref. 043CLN	• 6 BLACK CAPS NORMAL SIZE (ONLY FOR LABORATORY USE)	€ 11.00
	Ref. 040CRNDR8	• 6 LIGHT GREEN CAPS NORMAL SIZE (STRONG RETENTION) (undersized internal diameter for worn out sphere or for Ø 2.2 SPHERE)	€ 26.00
	Ref. 040CRNSNDR8	• 6 PINK CAPS NORMAL SIZE (SOFT RETENTION) (undersized internal diameter for worn out sphere or for Ø 2.2 SPHERE)	€ 26.00
	Ref. 049PCNDR8	• 6 ORANGE CAPS NORMAL SIZE (VERY ELASTIC RETENTION) (undersized internal diameter for worn out sphere or for Ø 2.2 SPHERE)	€ 26.00
	Ref. 042AMN	• 2 MONO OT BOX NORMAL SIZE (CASTABLE HOUSINGS FOR CAPS)	€ 19.00
	Ref. 041CAN B	• 2 STAINLESS STEEL HOUSINGS NORMAL SIZE	€ 18.00
	Ref. 041CTN B	• 2 TITANIUM HOUSINGS NORMAL SIZE	€ 27.00
	Ref. 101PPN	• 2 PLASTIC POSITIONER NORMAL SIZE	€ 5.00
	Ref. 044PPN	• 2 NORMAL TRANSFER PIVOTS	€ 22.00
	Ref. 100SAN	• 10 ALUMINIUM SPACERS FOR OT BOX BARS NORMAL SIZE	€ 7.00

MICRO Ø 1.8 mm:





	Ref. 055SCM	• 4 CASTABLE SPHERES MICRO SIZE	€ 19.00
	Ref. 151BAM	• 2 CASTABLE PRESECTIONED BEVELED BARS MICRO SIZE	€ 19.00
	Ref. 150BPM	• 2 CASTABLE PRESECTIONED BARS MICRO SIZ	€ 19.00
	Ref. 048CAM	• 6 SILVER CAPS MICRO SIZE (EXTRA-RESILIENT) ELASTIC AND GUMMY	€ 26.00
	Ref. 048COM	• 6 GOLD CAPS MICRO SIZE (EXTRA-RESILIENT) SLIGHTLY ELASTIC	€ 26.00
	Ref. 049PCM	• 6 GREEN CAPS MICRO SIZE (VERY ELASTIC RETENTION)	€ 26.00
	Ref. 060CRM AY	• 6 YELLOW CAPS MICRO SIZE (EXTRA SOFT RETENTION)	€ 26.00
	Ref. 040CRM SN	• 6 PINK CAPS MICRO SIZE (SOFT RETENTION)	€ 26.00
	Ref. 040CRM	• 6 WHITE CAPS MICRO SIZE (STANDARD RETENTION)	€ 26.00
	Ref. 040TCM	• 2 TITAN CAPS MICRO SIZE	€ 36.00
	Ref. 043CLM	• 6 BLACK CAPS MICRO SIZE (ONLY FOR LABORATORY USE)	€ 11.00
	Ref. 040CRMDR8	• 6 LIGHT GREEN CAPS MICRO SIZE (STRONG RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE)	€ 26.00
	Ref. 049PCMDR8	• 6 ORANGE CAPS MICRO SIZE (VERY ELASTIC RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE)	€ 26.00
	Ref. 042AMM	• 2 MONO OT BOX MICRO SIZE (CASTABLE HOUSINGS FOR CAPS)	€ 19.00
	Ref. 041CAM B	• 2 STAINLESS STEEL HOUSINGS MICRO SIZE	€ 18.00
	Ref. 041CTM B	• 2 TITANIUM HOUSINGS MICRO SIZE	€ 27.00
	Ref. 101PPM	• 2 PLASTIC POSITIONER MICRO SIZE	€ 5.00
	Ref. 044PPM	• 2 MICRO TRANSFER PIVOTS	€ 22.00
	Ref. 100SAM	• 10 ALUMINIUM SPACERS FOR OT BOX BARS MICRO SIZE	€ 7.00

NORMO-MICRO TOOLS:







	Ref. 085IAC	• INSERTION TOOL FOR CAPS - NORMAL AND MICRO SIZES (for metal housing and for inserting the caps in the mouth)	€ 22.00
	Ref. 74AC01	• KEY FOR PARALLELOMETER NORMAL SIZE	€ 22.00
	Ref. 74AC02	• KEY FOR PARALLELOMETER MICRO SIZE	€ 22.00
	Ref. 74GTN	• TECNO KEY FOR PARALLELOMETER NORMAL SIZE only for OT CAP TECNO (Gold Color)	€ 22.00
	Ref. 74GTM	• TECNO KEY FOR PARALLELOMETER MICRO SIZE only for OT CAP TECNO (Gold Color)	€ 22.00
	Ref. 091EC	• REMOVAL TOOL FOR RETENTIVE CAPS	€ 24,00
	Ref. 124ICP	• BLUE PLASTIC "MULTIUSE" INSERTION TOOL (for steel tool and for inserting the caps in the mouth)	€ 15,00

NEW

NORMAL-MICRO TOOLS continued:

	Ref. 080RCN	• REAMER TOOL FOR OT CAPS - NORMAL SIZE	€ 22.00
	Ref. 080RCM	• REAMER TOOL FOR OT CAPS - MICRO SIZE	€ 22.00
	Ref. 082ATN	• TOOL FOR TESTING CAP RETENTION NORMAL SIZE	€ 22.00
	Ref. 083ATM	• TOOL FOR TESTING CAP RETENTION MICRO SIZE	€ 22.00

NORMAL - MICRO ACCESSORIES AND 2.2mm:

	Ref. 155CC	• 6 CONNECTORS	€ 18.00
	Ref. 187CRS	• 2 CONCAVE SPHERES AVAILABLE: 1.8 - 2.2 - 2.5mm	€ 84.00
	Ref. 189SRS	• 2 SOLID SPHERES AVAILABLE: mm. 1,8	€ 80.00
	Ref. 100PD	• 10 PROTECTIVE DISKS	€ 7.00
	Ref. 038GCS	• 2 CASTABLE SLIDING BASES	€ 10.00
	Ref. 044PP22	• 2 MICRO TRANSFER PIVOTS 2.2 mm	€ 22.00

OT CAP ASSORTMENT KIT & CAP BOOKLET

KIT NORMAL

Ref. 074KAN

KIT MICRO

Ref. 074KAM

KIT CONTAINS:

- 8 WHITE CAPS - STANDARD RETENTION
- 8 PINK CAPS - SOFT RETENTION
- 6 YELLOW CAPS - EXTRA SOFT RETENTION
- 6 GREEN CAPS - VERY ELASTIC RETENTION
- 6 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC
- 6 SILVER CAPS - (EXTRA RESILIENT) ELASTIC AND GUMMY
- 4 LIGHT GREEN CAPS - STANDARD (FOR WORN OUT SPHERE)
- 4 ORANGE CAPS - ELASTIC (FOR WORN OUT SPHERE)
- 4 BLACK CAPS - LABORATORY USE ONLY
- 2 TITAN CAPS

€ 192.00

ASSORTMENT KIT NORMAL & MICRO

KIT CONTAINS:

- 4+4 WHITE CAPS - STANDARD RETENTION
- 4+4 PINK CAPS - SOFT RETENTION
- 2+2 YELLOW CAPS - EXTRA SOFT RETENTION
- 3+3 GREEN CAPS - VERY ELASTIC RETENTION
- 3+3 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC
- 3+3 SILVER CAPS - (EXTRA RESILIENT) ELASTIC AND GUMMY
- 2+2 LIGHT GREEN CAPS - STANDARD (FOR WORN OUT SPHERE)
- 2+2 ORANGE CAPS - ELASTIC (FOR WORN OUT SPHERE)
- 2+2 BLACK CAPS - LABORATORY USE ONLY
- 1+1 TITAN CAPS

€ 192.00

OT CAP ASSORTED RETENTIVE CAPS

NORMAL KIT

- 6 GREEN CAPS - VERY ELASTIC RETENTION
- 6 YELLOW CAPS - EXTRA SOFT RETENTION
- 6 PINK CAPS - SOFT RETENTION
- 6 WHITE CAPS - STANDARD RETENTION

€ 95.00

OT CAP ASSORTED RETENTIVE CAPS

MICRO KIT

- 6 GREEN CAPS - VERY ELASTIC RETENTION
- 6 YELLOW CAPS - EXTRA SOFT RETENTION
- 6 PINK CAPS - SOFT RETENTION
- 6 WHITE CAPS - STANDARD RETENTION

€ 95.00

Ref. 064ACN
NORMAL SIZE

Ref. 064ACM
MICRO SIZE



SPARE PARTS: OT STRATEGY


	Ref. 154PCS	• 4 OT STRATEGY STANDARD MALES	€ 19.00
	Ref. 054PLS	• 4 OT STRATEGY HIGH MALES	€ 19.00
	Ref. 154STY	• 2 CASTABLE STEADY	€ 9.00
	Ref. 047CSD AY	• 4 STRATEGY YELLOW CAPS FOR DUPLICATION TECHNIQUE (EXTRA SOFT RETENTION)	€ 26.00
	Ref. 047CSD SN	• 4 STRATEGY PINK CAPS FOR DUPLICATION TECHNIQUE (SOFT RETENTION)	€ 26.00
	Ref. 047CSD	• 4 STRATEGY WHITE CAPS FOR DUPLICATION TECHNIQUE (STANDARD RETENTION)	€ 26.00
	Ref. 047CLS	• 4 STRATEGY BLACK CAPS FOR DUPLICATION TECHNIQUE (FOR LABORATORY USE ONLY)	€ 11.00
	Ref. 047CSDDR8	• 4 AQUA STRATEGY CAPS FOR DUPLICATION TECHNIQUE (undersized internal diameter for worn sphere or 1.7mm sphere)	€ 26.00
	Ref. 045CSM AY	• 4 STRATEGY YELLOW CAPS FOR S/STEEL HOUSING (EXTRA SOFT RETENTION) (TO BE USED ONLY WITH HOUSINGS)	€ 21.00
	Ref. 045CSM SN	• 4 STRATEGY PINK CAPS FOR S/STEEL HOUSING (SOFT RETENTION) (TO BE USED ONLY WITH HOUSINGS)	€ 21.00
	Ref. 045CSM	• 4 STRATEGY WHITE CAPS FOR S/STEEL HOUSING (STANDARD RETENTION) (TO BE USED ONLY WITH HOUSINGS)	€ 21.00
	Ref. 045CLS	• 4 STRATEGY BLACK CAPS FOR S/STEEL HOUSING (FOR LABORATORY USE ONLY)	€ 11.00
	Ref. 045CSMDR8	• 4 AQUA STRATEGY CAPS FOR S/STEEL HOUSING (undersized internal diameter for worn sphere or 1.7mm sphere)	€ 21.00
	Ref. 046CSM	• 2 STRATEGY STAINLESS STEEL HOUSINGS	€ 20.00
	Ref. 101OTS	• 2 STRATEGY PLASTIC POSITIONER (For soldering stainless steel housings)	€ 5.00
	Ref. 047ACS	OT STRATEGY ASSORTMENT CAP KIT FOR DUPLICATION TECHNIQUE KIT CONTAINS: • 4 YELLOW CAPS - EXTRA SOFT RETENTION • 4 PINK CAPS - SOFT RETENTION • 4 WHITE CAPS - STANDARD RETENTION	€ 68.00
	Ref. 045ACS	OT STRATEGY ASSORTMENT CAP KIT FOR S/STEEL HOUSINGS KIT CONTAINS: • 4 YELLOW CAPS - EXTRA SOFT RETENTION • 4 PINK CAPS - SOFT RETENTION • 4 WHITE CAPS - STANDARD RETENTION	€ 55.00
	Ref. 044PPS	• 2 STRATEGY TRANSFER PIVOTS	€ 32.00
	Ref. 086ICS	• INSERTION TOOL FOR STRATEGY CAPS	€ 22.00
	Ref. 75AC04	• PARALLELOMETER KEY FOR STRATEGY CAPS	€ 22.00
	Ref. 081RCS	• REAMER TOOL FOR STRATEGY CAPS	€ 22.00

SPARE PARTS: OT BAR MULTIUSE



Ref. 022OBB	• 2 CASTABLE BARS OT BAR MULTIUSE	€ 16,00
Ref. 025CPB	• 4 CASTABLE BOXES OT BAR MULTIUSE	€ 18.00
Ref. 023CPA	• 4 POSITIONER CLIPS A	€ 16.00
Ref. 024CPB	• 4 POSITIONER CLIPS B	€ 16.00
Ref. 026CRR	• 4 PINK CLIPS (SOFT RETENTION)	€ 18.00
Ref. 027CRG	• 4 YELLOW CLIPS (MEDIUM RETENTION)	€ 18.00
Ref. 022RGO	• 2 CASTABLE CONNECTORS	€ 8.00
Ref. 029OIC	• TOOL FOR INSERTING OT BAR MULTIUSE CLIPS	€ 22.00
Ref. 028OCP	• PARALLELOMETER KEY FOR OT BAR MULTIUSE	€ 12.00

SPARE PARTS: OT VERTICAL










Ref. 072MCV	• 2 CASTABLE MALES • 2 PIN KEYS	€ 22.00
Ref. 072SCV	• 2 CASTABLE STEADY	€ 7.00
Ref. 072BCV	• 2 CERAMIC BARS	€ 6,00
Ref. 072CBV	• 4 WHITE CLIPS (STANDARD RETENTION)	€ 20.00
Ref. 072CVV	• 4 GREEN CLIPS (SOFT RETENTION)	€ 20.00
Ref. 072ICV	• 1 INSERTION TOOL FOR OT VERTICAL CLIPS	€ 22.00









SPARE PARTS: OT REVERSE 3



Ref. 120MR2	• 2 RETENTIVE MALES in TITANIUM + NYLON (Color WHITE/ STANDARD RETENTION)	€ 40.00
Ref. 041CAM A	• 2 S/STEEL HOUSINGS MICRO SIZE	€ 18.00
Ref. 121AOR	• 2 BRASS ANALOGS FOR ROOTS	€ 40.00
Ref. 121AAR	• 2 INOX IMPLANT ANALOGS	€ 50.00
Ref. 121MT	• 2 PLASTIC TRANSFER ANALOGS FOR IMPRESSION	€ 8.00

	Ref. 121AM	• 2 PLASTIC MANUAL TOOLS	€ 8.00
	Ref. 100PDR	• 10 PINK PROTECTIVE DISKS	€ 7.00
	Ref. 121PR	• 2 PLASTIC POSITIONERS FOR DUPLICATION	€ 8.00
	Ref. 124ICP	• PLASTIC BLUE INSERTING/HOLDER TOOL for caps	€ 15.00
	Ref. 124IMR	• METAL INSERTING TOOL FOR OT REVERSE/3 MALE	€ 15.00
	Ref. 72CBR	• S/STEEL PARALLELOMETER KEY FOR CASTABLE BAR	€ 30.00
	Ref. 72CFR	• BRASS PARALLELOMETER KEY FOR TITANIUM FEMALE	€ 30.00

SPARE PARTS: OT LOCK

	Ref. 882CS	• 1 SLIDING KEY (only for titanium - castable locking pin)	€ 14.00
	Ref. 882AF	• 1 FLEXIBLE RING	€ 18.00
	Ref. 882CB	• 1 LOCKING KEY	€ 20.00
	Ref. 882CPT	• 1 TITANIUM PREFABRICATED HOUSING (only for titanium - titanium adjustable locking pin)	€ 19.00
	Ref. 882PPC	• 1 CERAMIC PIN	€ 6.00
	Ref. 882FRC	• 1 CERAMIC SPACER AND CASTABLE RING (only for castable locking pin)	€ 6.00
	Ref. 882CLSCT	• 1 ADJUSTABLE SLIDING KEY (only for titanium adjustable locking pin)	€ 14.00
	Ref. 882DCC	• 9 CASTABLE SPACER RINGS (only for titanium adjustable locking pin)	€ 16.00

SPHERE REPAIR KITS:



AVAILABLE
1.8 - 2.2 - 2.5

Ref. 087CRS

CONCAVE REPAIR SPHERE

KIT CONTAINS:

- 2 CONCAVE SPHERES IN TITANIUM - TIN COATED
- 2 PINK CAPS SOFT RETENTION
- 1 INSERTING TOOL
- 1 GAUGE AND STRIP HOLDER

€ 100.00



Ref. 089SRS

SOLID REPAIR SPHERE

KIT CONTAINS:




- 2 SOLID SPHERES 1.8mm IN TITANIUM - TIN COATED
- 2 PINK CAPS SOFT RETENTION
- 2 PROTECTIVE DISKS
- 1 KEY

€ 100.00

TITANIUM ATTACHMENTS FOR OVERDENTURE OT REVERSE SYSTEM - OT KLASS:

OT REVERSE SPARE PARTS AND OT KLASS ARE IN STOCK
FOR DETAILED INFORMATION PLEASE E-MAIL US OR CALL OUR OFFICES

COMPLIMENTARY ITEMS

 <div>Height 14 cm</div>	Ref. 00P	• SMALL PARALLELOMETER DEVICE (FOR COURSES, ETC.)	€ 130.00
	◀ Base 12cm		
	Ref. 00BP	• MODEL HOLDER BASE (OPTIONAL)	€ 55.00
	Ref. 100AS	• 30 SILICON RINGS	€ 8.00
	Ref. NC	• Rhein'83 2 PART METAL BONDING COMPOSITE MATERIAL	€ 70.00


ADJUSTABLE DENTURE BOXES

IT IS POSSIBLE TO SET UP SIX DIFFERENT INTERIOR SIZES

	Ref. 0P6	KIT CONTAINS: • 6 ASSORTED BOXES COLOURS: BLUE - LIGHT BLUE - WHITE	€ 5.00

BROKEN SCREW EXTRACTOR KIT

BROKEN SCREW EXTRACTOR KITS ARE MANUFACTURED FOR MOST IMPLANT BRANDS.

	Ref. 780E FOR CORE VENT AND COMPATIBLE		
	Ref. 781E FOR BRANEMARK AND COMPATIBLE		
	Ref. 780	• 1 REVERSE CUTTING BUR	€ 50,00
	Ref. 780FL	• 1 CLAW REAMER BUR	€ 50,00

RHEIN'83 DIRECT SERVICES

Rhein'83 Direct Services

***A telephone service is available during office hours,
it furnishes direct answers for technical problems.
You can address to Rhein'83 for knowing the address of foreign distributors.***

Rhein'83 INTERNET: www.rhein83.com

E-mail: info@rhein83.com

Sales Dep.: angela.internationaltrade@rhein83.it

Technical support: gianni.storni@rhein83.it



Via E. Zago, 8
40128 (Bologna) Italy
Tel. (+39) 051 244510 - (+39) 051 244396
Fax (+39) 051 245238
Numero Verde 800-901172

<http://www.rhein83.com>
e-mail: info@rhein83.com



DISTRIBUTOR

