



 **SCHICK**
Dentalgeräte

T100

Made in Germany by Schick.

High – performance Turbine for precise finishing of all ceramic materials.



T100 – for professional long-term work

With the T100 you are in possession of a high – performance turbine matched to the special requirements of grinding all ceramic materials. It can be adapted to all customary milling units. The unique LED lighting (optional) facilitates the production of high-precision telescopes made of ceramic materials or even individual implant-abutments.

Your advantages Your profits

Enhanced tank capacity – about 1,3 l – minimizes re-filling with distilled water

You are saving precious time

A float gauge mounted in the tank is showing the filling quantity

Water does not spill over

Large pour hole for distilled water

Handy filling

The unit controls and regulates very precise pressure and water inflow to the turbine

Allows continuous running operation

Can be adapted to all customary milling units

An available milling unit can also be used in continuous operation to finish ceramic

Allows fitting with special SCHICK light head

Clear and absolute vision to the work object alleviates high-precision working

Nice design, robust, easy handling

Longevity



Technical Data

Length: 266 mm
Width: 155 mm
Height: 177 mm
Speed: 300.000 rpm
Operating air pressure: 2,5 bar
Weight: 3.2 kg

Scope of delivery

Turbine hand piece inclusive lubricant and tool to insert the burs
Control unit with integrated water tank about 1300 ml
Foot switch
All tubes with their connections
Measuring cup and funnel tube
Art.-N° 2640/1